Automated Impella Controller

**USER MANUAL
GEBRAUCHSANWEISUNG
MANUEL D'UTILISATION
MANUAL DEL USUARIO
MANUALE D'USO
ANVÄNDARHANDBOK
BRUGERVEJLEDNING
GEBRUIKERSHANDLEIDING
BRUKSANVISNING
INSTRUKCJA OBSŁUGI
UŽIVATELSKÝ MANUÁL**

ABIOMED
Recovering hearts. Saving lives.

TABLE OF CONTENTS

ENGLISH	3
DEUTSCH	34
FRANÇAIS.....	66
ESPAÑOL.....	98
ITALIANO	130
SVENSKA	162
DANSK	192
NEDERLANDS	224
NORSK	256
POLSKI.....	286
CESKY.....	318

ENGLISH

INTENDED AUDIENCE

This manual is intended for use by cardiologists, surgeons, nurses, perfusionists, and cardiac catheterization lab technicians who have been trained in the use of the Impella® Circulatory Support System.

WARNINGS

Use of the Impella System by trained and experienced practitioners has been associated with improved outcomes. Consequently, the first use of Impella® should be preceded by the completion of a contemporary Abiomed Impella training program and include on-site proctoring during the first use by Abiomed clinical support personnel certified in the use of Impella.

During defibrillation, do **NOT** touch the Impella Catheter, cables, or Automated Impella Controller.

Power the Automated Impella Controller using its internal battery if the integrity of the protective earth conductor is questionable.

Lithium-ion battery replacement by inadequately trained personnel could result in excessive temperatures, fire, or explosion. Only technicians authorized by Abiomed should remove or change the battery.

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

No modification of this equipment is allowed.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the electromagnetic compatibility (EMC) information provided in this manual.

During transport, the Automated Impella Controller may be exposed to stronger electromagnetic disturbance than during in-hospital use. Strong electromagnetic disturbance may cause the Automated Impella Controller to display soft button menu selections that were not selected by the user. Operators should be aware that, under these conditions, the operating parameters are not affected. No user intervention is required. Monitor Impella Catheter flow and patient hemodynamics to confirm normal operation. The condition will resolve itself once the Automated Impella Controller is no longer exposed to the disturbance.

Portable and mobile RF communications equipment can affect medical electrical equipment.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Use of cables, other than those sold by Abiomed, may result in increased emissions or decreased immunity of the Automated Impella Controller.

The Automated Impella Controller uses RFID (radio frequency identification) to identify and communicate with the purge cassette. Other equipment may interfere with the Automated Impella Controller even if that other equipment complies with CISPR emission requirements.

To prevent the risk of explosion, do **NOT** operate the Impella System near flammable anesthetics.

Do **NOT** use damaged or contaminated connector cables.

CAUTIONS

To prevent malfunction of the Automated Impella Controller, avoid long-term exposure to direct sunlight and excessive heat (40 °C).

To prevent overheating and improper operation, do **NOT** block the cooling vents of the Automated Impella Controller while it is operating.

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella Controller will operate for at least 60 minutes after the batteries have been fully charged.

Minimize exposure of Impella System components to sources of electromagnetic interference (EMI). Exposure to sources of EMI, such as cell phones and two-way radios, may cause operational interference. To clear interference, either increase the distance between system components and the EMI source or turn off the EMI source.

Operation of Impella System components may interfere with the operation of other devices. If interference occurs, increase the distance between the device and system components.

Have a backup Automated Impella Controller, purge cassette, connector cable, and Impella Catheter available in the unlikely event of a device failure.

Do **NOT** use the bed mount as a handle.

INTENDED USE (EU)

The Automated Impella Controller is intended for exclusive use of Impella® catheters and accessories. The controller may only be connected with devices which are listed in these instructions.

CONTRAINDICATIONS (EU)

There are no contraindications for use of the Automated Impella Controller, which operates the Impella pumps. Please see the respective catheter user manuals for indications, contraindications and possible complications.

OVERVIEW

REUSABLE SYSTEM COMPONENTS

The Impella System consists of the following reusable components:

- Automated Impella Controller (AIC)—The AIC provides an interface for monitoring and controlling the function of the Impella Catheters. The controller provides a fluid purge to the Impella Catheters. The controller provides backup power when the Impella System is operated away from AC power. The controller weighs 26 lbs (11.8 kg) and can operate on its internal battery for at least 60 minutes when fully charged.
- AIC cart—holds the controller. The cart has wheels for easy transport of the controller and a storage basket.

SINGLE-USE SYSTEM COMPONENTS

The Impella System also includes the following single-use components:

- Impella Catheter—The Impella Catheter is a microaxial blood pump that is operated by the AIC. There are different types of Impella Catheters, each of which has a programmable memory chip in the red (or blue) Impella plug. The AIC reads the memory chip to determine which type of Impella Catheter is connected and operate it as intended.


Single Use Components, continued

- **Purge cassette**—The purge cassette delivers rinsing fluid to the Impella Catheters. The purge fluid (typically 5% dextrose solution) flows from the purge cassette through the catheter to the microaxial blood pump to prevent blood from entering the motor.
- **Connector cable**—The white connector cable connects the Impella Catheter to the AIC. Clips on the cable are used to secure the purge tubing to the cable.

AUTOMATED IMPELLA CONTROLLER (AIC) FEATURES

IMPORTANT NOTE:

The underside of the AIC has a battery switch to turn on the batteries. This switch is turned off for shipping purposes. Before operating the AIC for the first time, make sure you turn this switch on. If the battery switch is not turned on, the AIC will not be able to operate on battery power.


Feature	Description
1 Display screen	Displays user information, including the labels for the soft buttons. (Display screen elements described in detail later in this section.)
2 Soft buttons	Display, open, and close menus. The function for each soft button is defined by labels adjacent to the button on the display screen; function changes depending on the screen. When the Impella Catheter is running, the default soft button labels are as follows: A) MUTE ALARM B) FLOW CONTROL C) DISPLAY D) PURGE MENU E) MENU
3 Power indicator	LED light above the selector knob; indicates the power status of the Automated Impella Controller. • Green light—controller is on and plugged into AC power or running on battery power • Amber light—controller is off but plugged into AC power • No light—controller is off and not plugged into AC power
4 Selector knob	Rotating push button; turn clockwise and counterclockwise to navigate through menu items; push to make a selection. Selector Knob Function <i>Rotate the selector knob on the controller to navigate through menu items. Push the selector knob to confirm your selection.</i>
5 Catheter plug	Connection point on the controller for the connector cable that connects to the Impella Catheter.
6 Purge disc	A flexible diaphragm on the purge cassette tubing used to monitor purge pressure and regulate purge flow.
7 Purge cassette	Contains the components for delivering the purge fluid; maintains the pressure barrier between the blood and the motor to prevent blood from entering the motor.
8 Purge cassette door	Spring-loaded door that opens to provide access to the purge cassette.

Feature	Description
1 Bed mount	Metal bracket on the back of the controller; attaches controller to the cart or bed
2 Purge cassette door release	Button located on the left side of the controller; press to open the purge cassette door
3 VGA OUT	Connection for connecting the controller to another monitor to slave the display
4 USB connector	Interface for data transfer by Abiomed maintenance or service personnel
5 Service	Connection used by Abiomed maintenance or service personnel
6 AC fuses	Electrical safety device in the event of current overload
7 AC plug	Connection point on the controller for the AC power cord
8 Power switch	Button that turns the controller on or off • ON: Press and hold the power switch for 2 seconds • OFF: (1) Disconnect the Impella Catheter from the Automated Impella Controller (2) Press and hold the power switch for 3 seconds (3) A pop-up confirmation box will appear (4) Press OK using the selector knob to confirm that the controller should be turned off NOTE: Do not hold down power switch for more than 30 seconds during operation--this initiates an emergency shutdown
9 Equipotential ground stud	Used to ground the Automated Impella Controller according to hospital procedures
10 Ethernet jack	Connection for downloading data

HOME SCREEN

The home screen displays operating parameters and information for the entire Impella® System.

Display Element	Description
1 Alarm window	The alarm window displays up to 3 alarms simultaneously, in order of priority from top to bottom. For each alarm, the alarm window displays: - Alarm header – displayed in the left column; window is color-coded red for critical alarms, yellow for serious alarms, white for advisory notifications, gray for resolved alarms - Detailed text – up to 3 lines of instructions for resolving the alarm condition are displayed in the right column of the alarm window next to the alarm header and subhead information
2 Catheter serial number	Displayed in the upper left of the display screen if a catheter is connected to the controller.
3 System date and time	The current date (YYYY-MM-DD) and time (24-hour format; HH:MM) are displayed in the upper center of the screen display. (In this example it is 2019 April 5 at 3:56 pm.)
4 AIC serial number and software version	The AIC serial number and current software version running on the console are displayed in the upper right hand corner of the display screen.
5 System power area	System power information is located at the top of the display screen, to the right. Battery status – bar within battery symbol indicates the overall remaining capacity of the batteries - Full green bar for fully charged battery - Partial green bar for battery that is at least 50% charged - Partial yellow bar for battery that is between 16% and 50% charged - Partial red bar for battery that is less than or equal to 15% charged - Moving gray bar for battery that is in charging mode - Numeric percentage of battery power remaining displayed to the left of the battery icon AC plug indicator AC plug indicator - White plug indicates that the controller is running on AC power - White plug with a red X indicates no AC power detected and the controller is running on battery power
6 Mute alarm indicator	Displayed in place of the words “MUTE ALARM” when an alarm is silenced. - Yellow bell with red X displayed when an alarm is muted - Not displayed when an alarm is active (but not muted) or when there are no active alarms

Display Element	Description
7 Soft button labels	The soft buttons on the Automated Impella Controller have corresponding labels adjacent to them on the display screen. These labels change depending on the type of screen displayed. MUTE ALARM • Mutes (silences) active alarms FLOW CONTROL (or NEXT) • FLOW CONTROL – For setting the current flow rate or performance level for the Impella Catheter • NEXT – Advances to the next screen DISPLAY (or BACK) • DISPLAY – Brings up the Display menu for viewing waveforms and navigating to other screen displays • BACK – Returns to the previous screen PURGE MENU (or EXIT) • PURGE MENU – Brings up the Purge System menu for changing the purge fluid, purge cassette, or purge system; de-airing the purge system; or transferring to the standard configuration • EXIT – Exits the current procedure MENU (or CANCEL) • MENU – Brings up a menu of options related to controller settings, alarm history, and starting a case • CANCEL – Exits out of current menu
8 Purge system area	Information about the purge system is displayed to the right of the flow area at the bottom of the display screen. Purge system marquee – scrolls from left to right when purge system is operating • Slow scrolling represents normal purge flow rate • Fast scrolling represents bolus flow rate and priming flow rate Y connector icon (for Impella 2.5 and Impella CP®) • Appears above the purge system marquee when the Impella System is configured using the Y connector in the set-up configuration Purge flow • Current purge flow displayed in mL/hr below the purge system marquee if the purge flow is known • Not displayed when the purge system is stabilizing, when there is no purge cassette, or when the procedure has not yet started Purge pressure • Current purge pressure (pressure of the purge fluid delivered through the catheter to the motor) displayed in mmHg below the purge flow
Note: Purge System Stabilization	<i>The purge system must stabilize after case start, a purge procedure, or resolution of a purge alarm. During this time, it may take up to 3 minutes for purge system information to display on the screen.</i>
9 Flow area	Information about Impella Catheter flow is displayed in the lower left corner of the display screen. Max / Min • Max / Min displays the range for the flow rate Current flow rate • Mean catheter flow displayed in liters per minute (L/min) – the numbers appear in <i>white</i> if the catheter position is correct; <i>yellow</i> if the catheter position is incorrect or unknown • If the system is unable to calculate flow, a yellow triangular caution icon is displayed with the message “Flow Calculation Disabled” Catheter operation icon • The circular catheter operation icon rotates when the Impella Catheter is running
10 Central display area	On the home screen, the central display area displays a heart pictogram and Impella Catheter position indicator message. Heart pictogram appears in the center of the home screen display. • Provides a visual representation of the current Impella Catheter position • Overlaid with a translucent yellow “?” when the controller cannot determine catheter position, position is wrong, or placement monitoring is suspended or disabled. Impella Catheter position indicator message displayed to the left of the heart icon. • Displays “Impella Position OK” in green when catheter position is correct • Displays “Impella Position Unknown” in yellow when catheter position is unknown • Displays “Impella Position in Ventricle” in yellow when catheter is in the ventricle • Displays “Impella Position Wrong” in yellow when catheter position is incorrect • Displays “Placement Monitoring Suspended” in yellow when there is a fault in the sensor • Displays “Placement Monitoring Disabled” in yellow when you turn off placement monitoring through the menu • Displays “Placement Monitoring Not Available” in yellow when a Impella CP with SmartAssist or Impella 5.5 with SmartAssist are used.

PLACEMENT SCREEN

The placement screen displays real-time operating data for the system. The screen displays the placement signal and motor current waveforms as well as the maximum/minimum and average values for each waveform in the central display area of the screen. Use the **DISPLAY** soft button to navigate to the placement screen.

Placement Signal Waveforms

The placement signal waveform displays a pressure measurement that is useful for determining the location of the open pressure area (Impella 2.5® and Impella CP®) or differential pressure sensor (Impella 5.0® and Impella LD®) of the catheter with respect to the aortic valve. The placement signal is used to verify whether the Impella Catheter is in the aorta or in the ventricle by evaluating the current pressure waveform as an aortic or ventricular waveform (Impella 2.5®, Impella CP®) or by evaluating the differential pressure as pulsatile or flattened (Impella 5.0® and Impella LD®). The scale for the placement signal waveform is displayed to the left of the waveform. The default scaling is 0–160 mmHg (Impella 2.5, Impella CP) or -20–60 mmHg (Impella 5.0 and Impella LD). For Impella 2.5 and Impella CP, it can be adjusted in 20 mmHg increments, with a minimum upper limit of 100 mmHg and a maximum upper limit of 240 mmHg. For Impella 5.0 and Impella LD, the maximum display range is -60–100 mmHg.

To the right of the waveforms is a display that labels the waveform, provides the units of measurement, and shows the maximum and minimum values and the average value from the samples received. At the bottom of that window is the time scale, which you can set by pressing the **DISPLAY** soft button.

Motor Current Waveform

Motor current is a measure of the energy intake of the Impella Catheter motor. The energy intake varies with motor speed and the pressure difference between the inlet and outlet areas of the cannula. Motor current provides information about the catheter position relative to the aortic valve. When the Impella Catheter is positioned correctly, with the inlet area in the ventricle and the outlet area in the aorta, the motor current is pulsatile because the pressure difference between the inlet and outlet areas changes with the cardiac cycle. When the inlet and outlet areas are on the same side of the aortic valve, the motor current will be dampened or flat because there is little or no pressure difference between the inlet and outlet areas.

The scale for the motor current waveform is displayed to the left of the waveform. The default scaling is 0–1000 mA. It is adjustable in 100 mA increments for the Impella Catheter, with a minimum difference between upper and lower limits of 200 mA and a maximum difference of 1000 mA.

To the right of the waveform is a display that labels the waveform, provides the units of measurement, and shows the maximum and minimum values and the average value from the samples received. You can set the time scale at the bottom of that window by pressing the **DISPLAY** soft button.

PURGE SCREEN

The purge screen displays purge system data. In the central display area of the screen, the purge flow rate and purge pressure are plotted as a function of time. To the right of the plots, the current purge flow rate and purge pressure are displayed. Use the **DISPLAY** soft button to navigate to the purge screen.

Purge Flow

The purge flow rate delivered by the purge cassette is displayed in mL/hr. The standard scale for the purge flow (0–30 mL/hr) is displayed to the left of the purge flow plot. The maximum value on this scale can be adjusted from 20 mL/hr to 200 mL/hr in increments of 10 mL/hr.

To the right of the plot is a display that labels the plot and shows the most recent value update. You can set the time scale at the bottom of the window by pressing the **DISPLAY** soft button.

A purge flow change notification can be enabled to indicate when the purge flow rate increases or decreases by 2.5 mL/h. The message is intended to aid patient management by alerting the clinician to changes in the rates of dextrose and heparin infusion through the purge fluid. The alarm clears when you press the **MUTE ALARM** button. This alarm is disabled by default. To enable this alarm, press **MENU**, select Settings/Service, and select Enable Purge Flow Change Notification.

Purge Pressure

The Automated Impella Controller regulates purge pressure, the pressure of the purge fluid delivered through the catheter to the motor. The purge pressure generated by the purge cassette is displayed in mmHg. The standard scale for the purge pressure (0–1500 mmHg) is displayed to the left of the purge pressure plot. The maximum value on this scale can be adjusted from 100 mmHg to 2000 mmHg in increments of 100 mmHg. The purge pressure in the system is set to an ideal pressure between 300–1100 mmHg and the purge flows are between 2–30 mL/hr. An alarm appears if purge pressure falls below 300 mmHg or exceeds 1100 mmHg.

To the right of the plot is a display that labels the plot and shows the most recent value update. You can set the time scale at the bottom of the window by pressing the **DISPLAY** soft button.

INFUSION HISTORY SCREEN

The infusion history screen displays the infusion volume as well as the amount of heparin and dextrose infused each hour. The current time period is displayed at the top of the list. Use the **DISPLAY** soft button to navigate to the infusion history screen.

MOBILE OPERATION

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella Controller will operate for at least 60 minutes after the batteries have been fully charged.

The Automated Impella Controller (AIC) can be operated on internal battery power when it is not connected to AC power.

1. Disconnect the AIC from AC power.
2. The controller beeps once every 5 minutes to alert you that it is running on battery power and a white advisory notification appears in the alarms area on the screen. The AC power icon displays a red X through it.
3. When the controller is connected back to AC power, the white advisory notification turns gray and the battery power icon turns green.

TRANSPORT WITHIN THE HOSPITAL

Patients supported with the Impella System may require transport within the hospital for various reasons. Transport can be safe and simple for patients supported with Impella Catheters.

Considerations for transport within the hospital:

- The AIC and Impella Catheters are designed to operate on battery power for at least 1 hour.
- Confirm that the battery capacity displayed on the controller is 100%.
- If transport time might be longer than 1 hour, bring an extension cord or confirm that you will be able to connect the controller to AC power once you arrive at your destination.
- When rolling the AIC cart across a threshold, firmly grasp the cart handle and pull it over the threshold.
- Pay close attention to all system components and connections when rolling the AIC cart over thresholds and through elevator doors.
- Do not stress the connector cable from the controller to the Impella Catheter.

TRANSPORT BETWEEN HOSPITALS

During transport, the Automated Impella Controller (AIC) may be exposed to stronger electromagnetic disturbance than during in-hospital use. Strong electromagnetic disturbance may cause the AIC to display soft button menu selections that were not selected by the user. Operators should be aware that, under these conditions, the operating parameters are not affected. No user intervention is required. Monitor Impella Catheter flow and patient hemodynamics to confirm normal operation. The condition will resolve itself once the AIC is no longer exposed to the disturbance.

Guidelines for Patient Transport

The AIC is qualified for safe use by healthcare professionals to facilitate the transport of patients supported by the Impella System from one medical facility to another. To use the system during transport you must understand and follow the instructions in this manual. The transport team should include person(s) fully trained in the use of the AIC and the Impella Catheter.

Patients are placed on the AIC for circulatory support for periods compliant with the intended use duration of the Impella Catheter. If, during this time, a patient requires additional resources and specialized teams at another location (e.g., a larger facility such as a transplant center), the patient may be transferred safely to such a location using the AIC. Maintaining optimal patient hemodynamic status and correct Impella Catheter position are two key factors in managing patients supported with the Impella System during transport. Steps should be taken to eliminate or minimize any aspect of the transport that might adversely affect these factors.

The AIC is designed to operate for 60 minutes on battery power. Transport teams should take this into consideration when planning the transport. If the total transport time is expected to include more than 60 minutes during which the system will be disconnected from AC power, arrangements should be made to use a vehicle with a built-in DC to AC power inverter.

Important Transport Considerations

1. Planning is critical to success. Abiomed representatives can help with planning for transport. They can be contacted 24 hours a day at the Abiomed 24hr hot-line (see back cover).
2. The Automated Impella Controller should be fully charged prior to transport. Keep the Automated Impella Controller connected to AC power (or an AC inverter) whenever possible.
3. Do not stress the connector cable from the controller to the Impella Catheter. Such tension could move the catheter out of correct position and compromise patient circulatory support.
4. Carefully monitor purge pressures during changes in altitude.
5. The controller should be positioned to allow easy access to the display screen and soft buttons to view alarms and make any necessary changes.
6. Maintain ACTs between 160 and 180 seconds or at the level recommended by the physician responsible for the patient.

PREPARING FOR USE IN FLIGHT—ENABLE FLIGHT MODE

Flight Mode is available in the software for air transport. When active, Flight Mode disables the purge cassette RFID transmitter so the controller cannot detect the purge cassette. If any purge system alarms are triggered during transport, the transport team should inform the managing hospital upon arrival.


When preparing to use the system in-flight, select > **Menu** > Settings/Service > **Enable Flight Mode**. When Flight Mode is active, a white alarm appears on screen (see Figure below). During air transport, a grounding cable is needed to connect the controller to ground. The grounding cable is provided by the hospital.

Once transport is complete, select > **Menu** > Settings/Service > **Disable Flight Mode**.

TO ENABLE FLIGHT MODE

Select: > **Menu** > Settings/Service > **Enable Flight Mode**.

A white alarm appears at the top of the screen:

You are instructed to connect the controller to a ground. Upon arrival, you are instructed to *disable* Flight Mode.

Note: the AIC cannot detect the purge cassette while in Flight mode. Therefore, there are special purge system alarms that direct you to inform the team receiving the patient. You are advised to have the team address operations when Flight mode is disabled.

TO DISABLE FLIGHT MODE

Select: > **Menu** > Settings/Service > **Disable Flight Mode**.

The updated software permits the user to switch off the Purge Cassette Detection (RF-ID) to enable safe air transport.

DEFIBRILLATION

During defibrillation, do NOT touch the Impella Catheter, cables, or Automated Impella Controller.

If emergency defibrillation is required during support with the Impella System, universal safety precautions should be followed. To avoid the risk of electrical shock, caregivers should use caution during defibrillation to ensure they do not touch the Impella Catheter, cables, or controller.

ECG INTERFERENCE

Operating the Automated Impella Controller may cause interference with electrocardiogram (ECG) signals. Please check the electrode pads and leads for good fixation and contact. If interference persists, activate the 50/100 Hz band-elimination filter or the 60/120 Hz band-elimination filter (also known as notch filter) on your ECG device. The filter frequency will be based on the AC power frequency for the country in which you are operating the equipment. If your ECG device does not have the appropriate filters, disconnect the Automated Impella Controller temporarily from AC power to obtain an undisturbed signal. Please observe the battery status while running the Automated Impella Controller on battery power.

LATEX

The AIC was not made with natural rubber latex or synthetic derivatives of natural rubber latex.

Battery Power Note

If the Automated Impella Controller is allowed to discharge completely and the system shuts down due to low battery, the controller will need to charge for an extended period of time before it will turn on again.

ALARMS OVERVIEW

The Automated Impella Controller monitors various functions to determine whether specific operational parameters are within expected limits. When a parameter goes outside of its specified limits, the Automated Impella Controller sounds an alarm message that can be viewed on the display screen on the front of the controller. The alarm tone indicates the severity of the alarm. The alarm message on the display screen is color-coded for severity and provides details on the cause of the alarm and how to resolve the alarm. After muting an alarm, if another alarm occurs it will be displayed but only heard if it is a higher priority alarm than the one that was muted.

ALARM LEVELS

Alarms are divided into three levels of severity:

- **Critical** (red)
- **Serious** (yellow)
- **Advisory** (white)

Category	Description	Audible Indicator*	Visual Indicator
Advisory	Notification	1 beep every 5 minutes	Alarm header on white background
Serious	Abnormal situation. Prompt action needed.	3 beeps every 15 seconds	Alarm header on yellow background
Critical	High priority. Immediate action needed.	10 beeps every 6.7 seconds	Alarm header on red background

* Sound pressure of audible alarm indicators is >67 dBA

Alarm Delay

For some AIC alarms, there is a short delay between the triggered event and the audible annunciation and visual display of the alarm—and the display of the warning on-screen. See Alarm Delay Information table, directly below.

Alarm Delay Information	
Impella Defective	8 second delay
Impella Position Wrong	11±5 second delay
Controller Error	12±3 second delay
Emergency Shutdown Imminent	15±1 second delay
Battery Failure	28±8 second delay
Controller Failure	38±8 second delay
Battery Comm. Failure	40±10 second delay
Purge System Blocked	75±45 second delay

Alarm Display

The alarm window is located in the upper left region of the display screen on the front of the AIC. Alarms are listed in order of priority, with the highest priority alarm at the top. Up to three alarms may be displayed at one time. The colored background behind the highest priority alarm will alternate between two shades of that color. The white panel displayed to the right of the alarm header contains instructions for resolving the alarm condition. The instructions should be followed in the order given.

Mute Alarm Function

Pressing the **MUTE ALARM** button on the upper right of the AIC display screen will silence the audible alarm indicator for 2 minutes (for red or yellow alarms) or 5 minutes (for white advisory alarms). When an alarm is silenced, the words **MUTE ALARM** next to the button are replaced by the mute alarm indicator, a crossed-out bell icon.

Alarm History Screen

The alarm history may be accessed through the **MENU**. This screen contains a log of the alarms that occurred during the case. The AIC maintains a long-term log that is saved after the AIC is powered down or after a power failure, and this information may be downloaded by Abiomed personnel.

Alarms That Resolve On Their Own

Alarm History			FLOW CONTROL ▶
171000	2019-08-21 15:17:08	Resolved: Flight Mode Enabled	P-8
171000	2019-08-21 15:17:04	Flight Mode Enabled	DISPLAY ▶
171000	2019-08-21 15:14:17	Resolved: Suction	
171000	2019-08-21 15:13:59	Resolved: Impella Position Wrong	
171000	2019-08-21 15:13:52	Suction	
171000	2019-08-21 15:10:57	Impella Position Wrong	
171000	2019-08-21 15:10:16	Resolved: Impella Stopped. Retrograde Flow	
171000	2019-08-21 15:10:13	Impella Stopped. Retrograde Flow	PURGE MENU ▶

The audible indicator will shut off if an alarm condition is resolved before you press **MUTE ALARM**. The visual message, however, will continue to be displayed, with the alarm header on a gray background, for 20 minutes or until you press **MUTE ALARM**. This allows you to identify the alarm that occurred.

ALARM MESSAGES-SUMMARY

Severity	Alarm Header	Action	Cause
Critical Alarms	Air in Purge System	The purge system has stopped. Press the Purge Menu soft key, then select De-Air Purge System.	There is air in the purge tubing.
	Battery Critically Low	Plug controller into AC power.	Battery power has 15% remaining capacity.
	Battery Failure	1. Plug controller into AC power. 2. Press switch located on the underside of the controller. 3. Switch to backup controller.	A battery switch is turned off or there is a malfunction of the switch.
	Battery Failure	Plug controller into AC power.	One of the batteries has failed.
	Battery Temperature High	Switch to backup controller.	Battery temperature is greater than 60°C.
	Complete Procedure	1. Follow the steps on the screen or 2. Exit the procedure	Complete Procedure serious alarm (yellow; see next page) is active and the user has not responded for an additional 2 minutes.
	Controller Failure	Switch to backup controller.	There is a problem with the controller electronics.
	Controller Failure	The purge system has stopped. Switch to backup controller.	The controller has detected a purge pressure sensor defect and has stopped the purge system.
	Emergency Shutdown Imminent	Release ON/OFF push button.	Power switch pressed for 15 seconds while Impella is still connected.
	Impella Disconnected	1. Check cable connection to console. 2. Check Impella connection to cable.	Running Impella Catheter disconnected.
	Impella Failure	Replace Impella.	There is a problem with the Impella Catheter motor.
	Impella Position In Ventricle (<i>Impella 2.5, Impella 5.0, Impella LD, and Impella CP</i>)	1. Reduce the P-level to P-2 until imaging is available. 2. Reposition with imaging. 3. Return to preferred P-level and confirm positioning with imaging.	Controller has detected that Impella Catheter is fully in the ventricle.
	Impella Position Wrong (<i>Impella 2.5, Impella CP, Impella 5.0, and Impella LD</i>)	1. Reduce the P-level to P-2 until imaging is available. 2. Reposition with imaging. 3. Return to preferred P-level and confirm positioning with imaging.	Controller has detected that Impella Catheter is in the wrong position.
	Impella Stopped	1. Restart Impella. 2. Replace Impella after 3rd unsuccessful restart attempt.	There may be a mechanical or electrical problem in the Impella Catheter.
	Impella Stopped	1. Replace white connector cable. 2. Switch to backup controller. 3. Replace Impella Catheter after 3rd unsuccessful restart attempt.	There may be a mechanical or electrical problem in the Impella Catheter.
	Impella Stopped Controller Failure	Attempt to restart catheter was unsuccessful. Switch to backup controller.	There is a problem with the controller electronics.
	Impella Stopped Retrograde Flow	To prevent retrograde flow, restart Impella or withdraw pump from ventricle.	Impella Catheter is not running possible retrograde flow through Impella Catheter.
	Impella Stopped Motor Current High	1. Restart Impella. 2. Replace Impella after 3rd unsuccessful restart attempt.	There is a problem with the Impella Catheter motor.
	Purge Disc Not Detected	Reinsert Purge Disc.	The controller is not detecting that the purge disc is clicked into the front of the controller.
	Purge Pressure High	1. Check purge system tubing for kinks. 2. Decrease concentration of dextrose in the purge solution.	Purge pressure is ≥ 1100 mmHg with the purge flow < 2 mL/hr.

Alarm Messages - continued

Severity	Alarm Header	Action	Cause
Critical Alarms	Purge Pressure Low	1. Check purge system tubing for leaks. 2. Increase concentration of dextrose in the purge solution. 3. Press the Purge Menu soft key then select Change Cassette and Bag.	Purge pressure has dropped below 300 mmHg with the purge flow ≥ 30 mL/hr for 30 seconds or longer.
	Purge Pressure Low (when Flight Mode enabled)	1. Check the purge system tubing for leaks. 2. Upon arrival at receiving hospital, notify managing team to address alarm condition once Flight Mode is disabled.	Purge pressure has dropped below 300 mmHg with the purge flow ≥ 30 mL/hr for 30 seconds or longer.
	Purge System Blocked	1. Check all purge system tubing for kinks or blockages. 2. Decrease concentration of dextrose in the purge solution.	Purge flow has dropped below 1 mL/hr. Kinked or blocked purge connecting tube. Kinked or blocked purge lumen in Impella Catheter.
	Purge System Failure	1. Replace Purge Cassette. Press the Purge Menu soft key, then select Change Cassette and Bag. 2. Switch to backup controller.	There is a problem with the purge cassette or purge unit driver.
	Purge System Failure (when Flight Mode enabled)	Upon arrival at receiving hospital, notify managing team to address alarm condition once Flight Mode is disabled.	There is a problem with the purge cassette or purge unit driver.
	Purge System Open	1. Check the purge system tubing for open connections or leaks. 2. Press the Purge Menu soft key then select Change Cassette and Bag.	Purge pressure has dropped below 100 mmHg for 20 seconds or longer.
	Purge System Open (when Flight Mode enabled)	1. Check the purge system tubing for open connections or leaks. 2. Upon arrival at receiving hospital, notify managing team to address alarm condition once Flight Mode is disabled.	Purge pressure has dropped below 100 mmHg for 20 seconds or longer.
	Retrograde Flow	Check for high afterload pressure.	Retrograde flow detected at high motor speed.
Serious Alarms	Battery Comm. Failure	Plug controller into AC power.	Loss of communication to the battery.
	Battery Level Low	Plug controller into AC power.	Battery has 50% remaining capacity.
	Battery Temperature High	1. Check controller for blocked air vents. 2. Switch to backup controller.	Battery temperature is greater than 50°C and less than or equal to 60°C.
	Complete Procedure	1. Follow the steps on the screen or 2. Exit the procedure	User has not responded to a de-air or purge procedure screen for more than 1 minute or a transfer to standard configuration screen for more than 5 minutes.
	Controller Error	Switch to backup controller.	There is a problem with the controller electronics.
	Defective Purge Cassette	Replace Purge Cassette. Press the Purge Menu soft key, then select Change Cassette and Bag.	There is a problem with the purge cassette hardware.
	Impella Catheter Not Supported	1. Replace Impella with supported catheter. 2. Contact Abiomed Service to upgrade Impella Controller.	The Impella Catheter is not supported to operate with the current version of controller software and/or hardware.
	Impella Defective	Do not use Impella. Replace Impella.	There is a problem with the Impella Catheter electronics.
	Impella Flow Low	1. Check for suction. 2. Check for high afterload pressure.	While Auto Mode is enabled, avg. flow is: • <1.0 L/min. for Impella 2.5 • <2.5 L/min for Impella CP or Impella CP with SmartAssist for more than 10 seconds.
	Impella Outflow Blocked	1. Confirm Impella position with imaging. 2. Pull Impella back 2 cm.	Flow to Impella Catheter outlet area obstructed.
	Impella Position Wrong	1. Confirm Impella position with imaging. 2. Pull Impella back 2 cm.	Controller has detected that the Impella Catheter is in the wrong position, with the outlet area too close to the aortic valve.
	Impella Stopped Controller Failure	Attempting to restart catheter. Locate backup controller.	Impella Catheter has stopped and controller is attempting to automatically restart the catheter.
	Placement Signal Not Reliable (Impella 2.5 and Impella CP)	Placement Monitoring and Suction are suspended. 1. Monitor patient hemodynamics. 2. Monitor Impella position with imaging. 3. Check the patient cable for kinks.	There is a problem with the Impella Catheter sensor signal.

Alarm Messages, continued

Severity	Alarm Header	Action	Cause
Advisory Alarms	Placement Signal Not Reliable (Impella 5.0 and Impella LD)	Placement and suction monitoring are suspended. Flow calculation is disabled. Refer to flow estimates in table below.	There is a problem with the Impella Catheter sensor signal.
	Placement Signal Not Reliable (Impella RP)	Monitor patient hemodynamics.	There is a problem with the Impella Catheter sensor signal.
	Purge Volume Critically Low	1. Press the Purge Menu soft key 2. then select Change Purge Fluid Bag.	There are 15 mL (in addition to 5% of the starting bag volume) or fewer remaining in the purge fluid bag.
	Reinstall Software	Software installation was unsuccessful. Reinstall software.	Software was not installed successfully.
	Suction	1. Reduce P-level. 2. Check Impella position. 3. Check filling and volume status.	Suction is detected.
	AC Power Disconnected	Controller is running on battery power.	AC power was disconnected.
	Audio Off	The audio for the following alarm has been disabled. <Alarm will be listed here>	User has disabled audio for Placement Signal Not Reliable, Purge Pressure High, Purge System Blocked, Suction, or Placement Signal Lumen Blocked alarm.
	Flight Mode Enabled	1. Connect controller to ground during air transport. 2. If equipped with Impella Connect, enable Flight Mode on module. 3. Upon arrival at receiving hospital, disable Flight Mode under MENU.	Flight mode has been enabled for transport.
	Impella Flow Reduced	1. Reduce P-Level 2. Check Impella position. 3. Check filling and volume status.	Motor speed has been reduced in response to suction.
	Impella Position Unknown (Impella 2.5 and Impella CP)	Impella position unknown due to low pulsatility, assess cardiac function.	Impella Catheter position unknown due to low pulsatility.
	Impella Position Unknown (Impella 2.5, Impella CP, Impella 5.0, and Impella LD)	Confirm Impella position with imaging.	Impella Catheter position unknown detected by algorithm.
	Placement Monitoring Disabled	Monitor Impella position with imaging or re-enable Placement Monitoring under MENU and Settings/Service.	User has disabled Placement Monitoring.
	Placement Signal Lumen Blocked (Impella 2.5 and Impella CP)	Placement and Suction Monitoring are suspended. Aspirate with syringe, then flush.	Placement signal pulsatility is low and placement signal amplitude is high (>130mmHg) for 30 consecutive seconds.
	Preventing Retrograde Flow	Impella P-level has increased to prevent retrograde flow. 1. Consider increasing target P-level. 2. For weaning, disable Retrograde Flow Control through MENU soft key.	Retrograde flow has been detected and minimum motor speed has been increased to more than target P-level.
	Purge Cassette Incompatible	Contact Abiomed Service to update Impella Controller.	Incompatible purge cassette RFID version.
	Purge Flow Decreased	The purge flow has decreased by 2.5 mL/hr or more. This is a notification only; no action is required.	Purge flow has decreased by ≥ 2.5 mL/hr.
	Purge Flow Increased	The purge flow has increased by 2.5 mL/hr or more. This is a notification only; no action is required.	Purge flow has increased by ≥ 2.5 mL/hr.

Alarm Messages, continued

Severity	Alarm Header	Action	Cause
	Purge Volume Low	1. Press the PURGE MENU soft key then select Change Purge Fluid Bag. 2. Follow the instructions to change the purge fluid.	There are 30 mL (in addition to 5% of the starting bag volume) or fewer remaining in the purge fluid bag.
	Surgical Mode Enabled	Impella pump stopped. Purge system running. 'Impella Stopped' alarm disabled. To exit this mode start Impella pump.	Surgical Mode has been enabled to silence "Impella Stopped" alarm at P-0.
	Transfer to Standard Configuration (Impella 2.5 and Impella CP)	Press the PURGE MENU soft key, then select Transfer to Standard Configuration.	Follow instructions or press MUTE ALARM to clear the alarm for 30 minutes.
	Optical Sensor Not Supported (Impella CP with SmartAssist and Impella 5.5 with SmartAssist)	Placement and Suction Monitoring are Not Available. 1. Monitor Impella position with imaging. 2. Switch to Controller which supports optical pressure sensor.	The optical sensor is not supported with the current version of the controller hardware. The pump will still operate without placement or suction monitoring.
	Unexpected Controller Shutdown	Switch to backup controller if condition persists.	Unexpected restart of controller due to software of hardware failures.

PURGE CASSETTE PROCEDURES

There are three procedures for maintaining the Impella Catheter purge system:

- Change Cassette & Bag
- Change Purge Fluid Bag
- De-Air Purge System


Each procedure can be accessed using the **PURGE MENU** soft button. The purge cassette procedures are discussed below.

CHANGE CASSETTE & BAG

The Purge Cassette may need to be replaced for extended use of the Impella system.

Perform the following steps to change both the Purge Cassette and Purge Fluid:

1. Press **PURGE MENU** and select "Change Cassette & Bag" from the menu.
2. Select **START** to begin the cassette and fluid change process.

3. When prompted by the controller, disconnect the luer from the Impella catheter.

Note on Connecting Purge Tubing to the Catheter:

If included, remove and discard the y-connector and connect the yellow luer on the purge tubing directly to the yellow check valve on the Impella Catheter.

4. Open the Purge Cassette door by pressing the button on the left side of the console. Remove and discard the old cassette and purge fluid bag.
5. Open the new Purge Cassette. Spike the new purge fluid bag with the new purge cassette tubing. Select **NEXT** to continue.
6. Insert the new purge cassette into the controller. Be sure to slide the purge disc into place and extend the purge tubing through the gap in the purge cassette door when you close the door.
7. Confirm the luer is disconnected. Press **NEXT** to proceed to prime the purge cassette.
8. Update the purge fluid information.
 - a. To select the default purge fluid values displayed on the screen, select **CONFIRM**.
 - b. To change the purge fluid information, select **EDIT**. Then use the soft keys to navigate selections and edit values. Select **DONE** to complete editing.
9. When Steps 1 through 8 are complete, connect the luer from the new purge cassette to the Impella catheter.

CHANGE PURGE FLUID BAG

Note on Purge Solution Bottles:

If the purge solution is supplied in bottles, open the vent on the purge fluid spike and follow the same procedure as if supplied in bags.

Perform the following steps to change only the purge fluid.

1. Press **PURGE MENU** and select "Change Purge Fluid Bag."
2. Select **START** to begin the purge fluid change process.
3. When prompted by the controller, remove the old purge bag and replace by spiking the new purge fluid bag. Select **NEXT** to advance to the next step.
4. Update the purge fluid information.
 - a. To select the default purge fluid values displayed on the screen, select **CONFIRM**.
 - b. To change the purge fluid information, select **EDIT**. Then use the soft keys to navigate selections and edit values. Select **DONE** to complete editing.

- When prompted by the controller, disconnect the luer from the Impella Catheter. The controller will automatically prime the tubing, which will flush the fluid from the last bag out of the purge cassette tubing.

a. To skip the flush select **SKIP PRIME**.

Note:

The instructions to disconnect the luer and to automatically prime the tubing only occurs if the user changed the purge fluid dextrose or heparin concentrations.

- When prompted by the controller, connect the yellow luer from the purge cassette to the Impella catheter.

Note on Priming Purge Cassette and Disconnecting Luers:

The system only prompts you to remove luers to prime and flush tubing if the heparin or dextrose concentrations were changed or air was detected. You can only select **SKIP PRIME** if no air was detected.

DE-AIR PURGE SYSTEM

Perform the following steps to de-air the purge system:

- Press **PURGE MENU** and select "De-Air Purge System."


Note on De-Air Procedure:

You may run the de-air procedure after changing the dextrose concentration to decrease the amount of time it takes for a change to occur.

- Select **START** to begin the de-air process.
- Make sure that the purge fluid bag is NOT empty or inverted and that the tubing is NOT kinked. Select **NEXT** to continue.
- Disconnect the purge tubing from the Impella Catheter.
- Confirm that no air remains in the purge tubing. If air remains, press **BACK** to repeat the air removal process.
- Connect the purge tubing to the luer on the Impella Catheter to complete the de-air procedure.

AIR DETECTED ALERT

During any of the purge system processes above, the controller automatically monitors for air in the system. If air is detected in the system, the controller provides an alert to disconnect the luer as shown below. Once the luer is disconnected, the controller automatically de-airs the purge system.

TROUBLESHOOTING THE PURGE SYSTEM

PURGE PRESSURE LOW

Note on Purge Pressure

The optimal purge pressure is different for every Impella Catheter. Purge pressure can range from 300 mmHg to 1100 mmHg. While purge pressure varies during operation, the Automated Impella Controller automatically maintains purge pressure within an acceptable range for each Impella Catheter.

Perform the following steps to troubleshoot the Purge Systems

- Inspect the purge system for leaks.
- If there are no leaks, change to a purge fluid with a higher dextrose concentration. To do this, open the **PURGE MENU** and select "Change Purge Fluid Bag." Follow the instructions on the screen. (Refer to "Purge Cassette Procedures" earlier in this section of the manual.)
- If the pressure stabilizes, no other action is required.
If the purge pressure is not stable, proceed to Step 4.
- If the low purge pressure alarm remains unresolved for more than 20 minutes, there may be a problem with the purge cassette. Replace the purge cassette. (Refer to the "Change Cassette & Bag" instructions, earlier in this section.)

PURGE PRESSURE HIGH AND PURGE SYSTEM BLOCKED

If the purge pressure exceeds 1100 mmHg, the Automated Impella Controller displays the "Purge Pressure High" alarm. If the purge flow stops completely, the controller displays the "Purge System Blocked" alarm.

For either event, follow these steps:

- Inspect the purge system and check the Impella Catheter for kinks in the tubing.
- Check the glucose concentration of the purge fluid. Decrease the concentration to 5% if current concentration is higher
- Replace the purge cassette using the Change Cassette and Bag procedure shown earlier in this section.
- Monitor Motor Current.

Note on Unresolved Purge Pressure High Alarm:

If not resolved by the recommendations provided, high purge pressure—which triggers the "Purge Pressure High" alarm message—could be an indication of a kink in the Impella Catheter. In this case, the motor is no longer being purged and may eventually stop. Clinicians should monitor motor current and consider replacing the Impella Catheter whenever a rise in motor current is seen.

PURGE SYSTEM OPEN

Note on Purge System Open Alarm

This alarm may occur if purge pressure is less than 100 mmHg for 20 seconds or longer.

Perform the following steps to check the Purge System:

- Inspect the purge system for leaks.
- If no leaks are visible, there may be a problem with the purge cassette. Replace the purge cassette. (Refer to the "Change Cassette & Bag" instructions, earlier in this section.)

SYMBOLS

	Caution: Read, understand, and follow instructions for use
	Defibrillator-proof type CF equipment
	Keep dry
	Storage temperature (eg, 10°C to 30°C)
	Declares conformity with Directive 93/42/EEC for medical devices or with Regulation (EU) 2017/745 on medical devices, and with Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment
2014-10-01	Date of manufacture (eg, October 1, 2014)
	Protect from sunlight
	Symbol for lot designation; the manufacturer's lot designation must be stated after the LOT symbol
REF 123456	Abiomed part number (eg, part number 123456)
SN 123456	Manufacturer's serial number (eg, serial number 123456)
Non Sterile!	The product is not sterile
2016-06-01	Use-by date (eg, use before June 1, 2016)
	Do not reuse

	Sterilized using ethylene oxide
	Medical Device
	Electric scrap; must be disposed of separately. Must not be disposed of as domestic waste.
	Protective Earth
	ON / OFF
	Alternating current (AC) only
	Equipotentiality
	Fuse
	Non-ionizing electromagnetic radiation
	USB port
	CAT 5 Port (Ethernet)
	MR Unsafe
Do Not Flush	Do NOT flush
Glucose	Use glucose in the purge fluid
NaCl	Use saline in the pressure bag; squeeze at the green arrows to flush

AIC MECHANICAL SPECIFICATIONS

Parameter	Specification
Model Number	0042-0000-EU / 0042-0000-CA
Temperature	Operating: 10 °C to 40 °C (50 °F to 104 °F) Storage: -15 °C to 50 °C (5 °F to 122 °F)
Relative Humidity	Operating: 95% Storage: 95%
Atmospheric Pressure	Operating: 8000 ft (750 hPa) to -1000 ft (1050 hPa) Storage: 18,000 ft (500 hPa) to -1000 ft (1050 hPa)
Dimensions	Height: 351 mm (13.8 in) Width: 443 mm (17.4 in) Depth: 236 mm (9.3 in)
Dimensions – Packaged	Height: 508 mm (20.0 in) Width: 559 mm (22.0 in) Depth: 406 mm (15.0 in)
Weight	Maximum: 12.2 kg (26.8 lbs)
Weight Packaged	Maximum: 14.3 kg (31.5 lbs)
Maintenance and repair interval	12 months (Work must be performed by technicians authorized by Abiomed and have completed Abiomed's Service Training Certification Program.)

AIC ELECTRICAL SPECIFICATIONS

AC operation	100-240 V AC, 2A, 50-60 Hz
Internal battery operation	14.4 V DC (nominal); lithium ion
Characteristic values	Max. power consumption under load: 200 VA 9.7 fuses: 2 Amp. 250 V. 5 mm x 20 mm, slow-blow fuses Running time without AC power, with fully charged batteries: At least 60 minutes (charging duration of at least 5 hours)
Electrical system	Installation in accordance with pertinent regulations is required for use in medical facilities (eg, VDE 0100, VDE 0107, or IEC stipulations). Observe country-specific regulations and national deviations.

NOTE Circuit diagrams are available upon request.

EQUIPMENT CLASSIFICATIONS

Type of protection against electric shock	IEC 60601-1: Class I degree of protection: CF defibrillation-proof and internally powered. Relies not only on basic insulation against shock but also includes additional protection. Accomplished by providing means for connecting the equipment to the protective earth conductor of the fixed wiring of the installation in a way that prevents accessible metal parts from becoming live if basic insulation fails.
Degree of protection against electric shock for Automated Impella Controller	Class I Equipment
Mode of operation	Continuous
Degree of protection against explosion hazard	Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Also not suitable for use in an oxygen-enriched atmosphere.
Degree of protection against harmful ingress of water	IEC 60529: IPX1 protected against dripping water.

EQUIPMENT DESIGN

The Automated Impella Controller conforms to the applicable requirements of the following standards:

- IEC 60601-1: 2012 Edition 3.1 Medical Electrical Equipment Part 1: General Requirements for Basic Safety and Essential Performance
- CSA C22.2#60601-1 (2014) Ed:3 Medical Electrical Equipment Part 1: General Requirements for Basic Safety and Essential Performance
- AAMI ES60601-1:2005 +C1:A2 Medical Electrical Equipment Part 1: General Requirements for Basic Safety and Essential Performance
- IEC 60601-1-2:2014 Edition 4, Medical Electrical Equipment - Part 1-2: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirements and Tests
- IEC 60601-1-6:2010, AMD1:2013 Medical Electrical Equipment - Part 1-6: General Requirements for Safety - Collateral Standard: Usability
- IEC 60601-1-8:2006, AM1:2012 Medical Electrical Equipment - Part 1-8: General Requirements for Safety - Collateral Standard: General Requirements, Tests and Guidance for Alarm Systems in Medical Electrical Equipment and Medical Electrical Systems
- IEC 62304:2015 Medical Device Software - Software Life-cycle Processes
- RTCA DO160G Environmental Conditions and Test Procedures for Airborne Equipment
- AIM 7351731 Medical Electrical Equipment and System Electromagnetic Immunity Test for Exposure to Radio Frequency Identification Readers

Complete list of standards available upon request.

FEDERAL COMMUNICATIONS COMMISSION NOTICE

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by Abiomed, Inc. could void the user's authority to operate this device.

ELECTROMAGNETIC COMPATIBILITY

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the electromagnetic compatibility (EMC) information provided in this document.

Portable and mobile RF communications equipment can affect medical electrical equipment.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Use of cables, other than those sold by Abiomed, may result in increased emissions or decreased immunity of the Automated Impella Controller.

The Automated Impella Controller uses RFID (radio frequency identification) to identify and communicate with the purge cassette. Other equipment may interfere with the Automated Impella Controller even if that other equipment complies with CISPR emission requirements.

The Automated Impella Controller (AIC) performs as intended when exposed to radiofrequency (RF) disturbances below 20 V/m. During transport, the AIC may be exposed to RF disturbances above 20 V/m, which could cause minor problems, such as intermittent displays of soft button menu selections, which have no effect on the operating parameters of the Impella support system, and will resolve readily once the disturbance ends. It could also potentially result in loss of support. Patients must be closely monitored at all times during transport.

Do not transport an Impella patient via commercial aircraft. Loss of support may occur aboard a commercial aircraft due to exposure to radiofrequency (RF) disturbances above the compliance level (<20 V/m) of the Automated Impella Controller.

NOTE: The EMC tables and other guidelines included in this manual provide information to the customer or user that is essential in determining the suitability of the equipment, or system, for the electromagnetic environment of use, and in managing the electromagnetic environment of use permit the equipment or system to perform to its intended use without disturbing other equipment and systems or non-medical electrical equipment. For the electromagnetic testing (detailed in the following tables), the AIC Essential Performance was specified as: during the entire testing period, the AIC continues to provide support to the patient.

Table 1. (Table 201) Guidance and Manufacturer's Declaration – Emissions, All Equipment and Systems

The Automated Impella Controller is intended for use in the electromagnetic environment specified below. The customer or user of the Automated Impella Controller should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Enforcement – Guidance
RF Emissions CISPR 11	Group 1 Class A	The Automated Impella Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonics IEC 61000-3-2	Class A	The Automated Impella Controller is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Flicker IEC 61000-3-3	Complies	
RTCA DO-16G Section 21.4, conducted emissions	Category M	
RTCA DO-16G Section 21.5, radiated emissions	Category B	

NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Table 2. (Table 202) Guidance and Manufacturer's Declaration – Immunity

The Automated Impella Controller is intended for use in the electromagnetic environment specified below. The customer or user of the Automated Impella Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	The relative humidity should be at least 5%.
Electrical Fast Transient/burst IEC 61000-4-4	±2 kV Mains ±1 kV for input/output lines	±2 kV Mains ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV Differential ±2 kV Common	±1 kV Differential ±2 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	> 95% dip for 0.5 cycle 60% dip for 5 cycles 30% dip for 25 cycles > 95% dip for 5 seconds	> 95% dip for 0.5 cycle 60% dip for 5 cycles 30% dip for 25 cycles > 95% dip for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Automated Impella Controller requires continued operation during power mains interruptions, it is recommended that the Automated Impella Controller be powered from an uninterruptible power supply or battery.
Power Frequency 50/60 Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.

Table 3. (Table 203) Guidance and Manufacturer's Declaration – Emissions Equipment and Systems that are Life-Supporting

The Automated Impella Controller is intended for use in the electromagnetic environment specified below.
The customer or user of the Automated Impella Controller should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be separated from the Automated Impella Controller by no less than the recommended separation distances calculated/listed below:
Conducted RF IEC 61000-4-6	10 Vrms 150 kHz to 80 MHz	10 Vrms	$d = 0.35\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	20 V/m	$d = 0.6\sqrt{P}$ 80 MHz to 800 MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and d is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey ^(a) , should be less than the compliance level in each frequency range. ^(b) Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Automated Impella Controller is used exceeds the applicable RF compliance level above, the Automated Impella Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Automated Impella Controller.

^(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Table 4. (Table 205) Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the AIC, Equipment and Systems that are Life-Supporting

The Automated Impella Controller is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Automated Impella Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the Automated Impella Controller as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Output Power of Transmitter (Watts)	Recommended Separation Distances for the Automated Impella Controller (m)		
	150 KHz to 80 MHz $d = 0.35\sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.5 GHz $d = 1.2\sqrt{P}$
0.01	0.04	0.06	0.12
0.1	0.11	0.19	0.38
1	0.35	0.6	1.2
10	1.11	1.9	3.8
100	3.5	6.0	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Immunity Test	Compliance Level
Avionics: RTCA DO-160G	
Conducted RF: Section 20.4	Category R: 10 MHz - 400MHz
Radiated RF: Section 20.5	Category T (c): 100 MHz - 8 GHz
(c) the AIC will not maintain its essential performance when subjected to Category R Levels (radiated RF at a field strength of 150 V/m).	

Immunity Test		
RFID AIM 7351731:2017		
RFID Specification	Frequency	Test Level (RMS)
ISO 14223	134.2 kHz	65 A/m
ISO/IEC 14443-3 (Type A)	13.56 MHz	7.5 A/m
ISO/IEC 14443-4 (Type B)	13.56 MHz	7.5 A/m
ISO/IEC 15693 (ISO 18000-3 Mode 1)	13.56 MHz	12 A/m
ISO/IEC 18000-7	433 MHz	3 V/m
ISO/IEC 18000-63 Type Ca	860-960 MHz	54 V/m
ISO/IEC 18000-4 Mode 1	2.45 GHz	54 V/m

RFID Transmitter / Receiver Specifications

Frequency	13.56 MHz
Receiver bandwidth	14 kHz
Effective radiated power	30 nW
Modulation	ASK

WHITE CONNECTOR CABLE

Length	2.5 m
Service life	Single use only

SLAVE MONITOR CONNECTION

If equipped with a VGA connector, the Automated Impella Controller can be connected to a monitor to display information from the controller, at a resolution of 800 x 600 pixels. The connection between the controller and the monitor can be made using a cable up to 20 feet in length.

CLEANING

- Clean the Automated Impella Controller keypad and display with either 70% isopropyl alcohol or soap and water. (NOTE: Be aware that soft buttons may be activated when you spray or wipe the display)
- Clean the Automated Impella Controller housing with mild detergent
- Sani-Cloth HB wipes may be used to clean the AIC.
- Do NOT clean with or expose any part of the clear sidearm of the Impella Catheter (e.g., infusion filter, pressure reservoir) to alcohol. Alcohol has been shown to cause cracks and leaks in these components.
Carefully read labels on common skin preps and lotions to avoid using any alcohol-containing products in the area of the infusion filter or pressure reservoir.
- Do NOT allow any fluids to enter the connector sockets.
- Clean the connector cable with 70% isopropyl alcohol.

Alcohol Warning

Do NOT clean the Impella Catheter infusion filter or pressure reservoir with alcohol and AVOID exposing these components to products containing alcohol.

STORING THE AUTOMATED IMPELLA CONTROLLER

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella Controller will operate for at least 60 minutes after the batteries have been fully charged.

Place the AIC on a horizontal surface to prevent falling.

- Connect the AC power cord to an AC outlet.
- The battery may be destroyed if the Automated Impella Controller is stored with a depleted battery.
- To keep the Automated Impella Controller battery charged, the controller should be plugged into an AC outlet. When plugged into an AC outlet, the controller battery will charge whether the controller is on or off.

Battery Power Note

If the Automated Impella Controller is allowed to discharge completely and the system shuts down due to low battery, the controller will need to charge for an extended period of time before it will turn on again.


IMPELLA SYSTEM LIMITED SERVICE WARRANTY

Impella System Limited Service Warranty (Europe)

Please contact your EU representative for warranty information.

INSPECTIONS, MAINTENANCE AND REPAIR

The Automated Impella Controller (AIC) is subject to service in accordance with §7 MPBetrV (German regulation). Service includes in particular maintenance and repair. Maintenance activities include the yearly inspections and maintenance required to ensure the continued safe and correct operation of the AIC. The corresponding work may only be carried out by technicians authorized by Abiomed.

A sticker on the device indicates the date for the next required inspection. This example shows a sticker indicating that inspection is required in June 2020.


The following inspections are required for the Automated Impella Controller:

- Inspection of labeling and instructions for use
- Visual inspection of the device and its accessories for any signs of damage
- Testing for electrical safety as per DIN VDE 751 or DIN EN 60 601
- Leakage current test
- Dielectric strength test
- Functional testing of all switches, keys, rotary knobs, sockets, and control lights on the device
- Checking battery operation

If defects become apparent during the technical safety inspections that could endanger patients, employees, or third parties, then the device must not be operated until the defects have been remedied by proper technical servicing.

PATIENT ENVIRONMENT

The Automated Impella Controller and the components of the Impella System are approved for use within the patient environment, as defined in IEC 60601-1-1 / IEC 60601-1: 3rd edition, and in the following figure.

PATIENT MANAGEMENT OVERVIEW

The following section addresses left-heart support with Impella 2.5, Impella CP, Impella CP with SmartAssist, Impella LD, Impella 5.0, or Impella 5.5 with SmartAssist heart pumps only. Please refer to the Impella RP User Manual for right-heart support. The information and instructions in this section are not intended to supersede established medical procedures concerning patient care. Best practices—as determined by the medical community—should always be observed. In each case, the clinician must determine whether the application of information provided is appropriate for the particular clinical setting.

PATIENT SELECTION FOR LEFT-VENTRICULAR SUPPORT

INITIAL PATIENT DECISION: PCI VS. CABG

Per the ESC/EACTS guidelines, there should be an initial overall decision around patient selection relative to the index procedure being CABG or PCI. The PCI indication is a medical decision and should be made by a heart team according to institution standards, current practice of medicine, and societal ESC/EACTS guidelines. This involves a multidisciplinary approach composed of clinical or non-invasive cardiologists, cardiac surgeons, and interventional cardiologists, as well as anaesthetists and other specialists, if deemed necessary.

Impella Patient Decision

For the use of the Impella 2.5 and the Impella CP, the patient would be deemed an appropriate candidate for high-risk PCI as defined by the inclusion/exclusion criteria contained in the heart pump's IFU.

GENERAL PATIENT CARE CONSIDERATIONS

- Use knee immobilizer as needed to maintain access site straight.
- Access site management should be done in accordance with hospital protocol, using aseptic technique.
- Assess access site for bleeding and hematoma.
- Monitor pedal pulses.
- To prevent the purge tubing from kinking, do not allow the red Impella plug to hang freely from the catheter and do not bend the catheter near the red Impella plug.
- Consider attaching the red Impella plug and catheter to a short armboard to prevent the catheter from kinking near the plug.
- When transferring a patient with the device in place: Be careful not to pull on the Impella Catheter when transferring a patient from one bed to another.
 - Do not raise the head of the bed to higher than a 30-degree angle.
 - Use care when moving or turning a patient; the Impella Catheter may move out of position and cause a positioning alarm.
 - ACT 160-180 seconds

RIGHT HEART FAILURE

Caregivers should monitor patients being supported by the left-ventricular Impella Catheter for signs of right heart failure:

- Reduced output from the Impella Catheter
- Suction alarms
- Elevated filling pressures (CVP)
- Signs of liver failure
- Elevated pulmonary pressures

If the patient is exhibiting signs of right heart failure, the clinical team should assess the need for a more durable form of support.

USE OF ECHOCARDIOGRAPHY FOR POSITIONING IMPELLA CATHETER IN THE LEFT VENTRICLE

BACKGROUND

Echocardiography is a commonly used tool for evaluating the position of the Impella Catheter relative to the aortic valve and other intraventricular structures post-placement. The best echocardiographic views for positioning the Impella Catheter in the left ventricle are a long axis transesophageal echocardiogram (TEE) or a parasternal long axis transthoracic echocardiogram (TTE). These long axis views allow you to see both the aortic valve and Impella Catheter inlet area.

Evaluate the position of the Impella Catheter if the Automated Impella Controller displays position alarms or if you observe lower than expected flows or signs of hemolysis. If the catheter does not appear to be correctly positioned, initiate steps to reposition it. The following illustrations identify the structures you would expect to see in transesophageal echocardiography (top) and transthoracic echocardiography (bottom).

In these illustrations, the Impella Catheter is positioned correctly; however, these depictions are stylized. In actual echocardiograms, the pigtail and inlet and outlet areas may not be seen as distinctly.

The graphics in this section depict the Impella 2.5 Catheter, but are representative of positioning for the Impella CP and Impella 5.0 as well. Positioning for the Impella LD and Impella 5.5 Catheters would also look similar, but without the pigtail.


Figure 1. Transesophageal Echocardiogram (TEE) of Impella Catheter


Figure 2. Transthoracic Echocardiogram (TTE) of Impella Catheter

FOUR IMPELLA CATHETER POSITIONS

There are four Impella Catheter positions you are likely to encounter when examining echocardiograms from patients supported with the Impella Catheter include:

- 1-Correct Impella Catheter position
- 2-Impella Catheter too far into the left ventricle
- 3-Impella Catheter inlet in the aorta
- 4-Impella Catheter in papillary muscle

The following pages describe each situation. Figures 3–6 illustrate transesophageal echocardiograms (TEE). Figures 7–10 illustrate transthoracic echocardiograms (TTE) of each situation.

1—CORRECT IMPELLA CATHETER POSITION

For optimal positioning of the Impella Catheter, the inlet area of the catheter should be 3.5 cm below the aortic valve annulus and well away from papillary muscle and subannular structures. The outlet area should be well above the aortic valve. If the Impella Catheter is correctly positioned, echocardiography will likely show the following, as depicted in Figure 1 and 3 (TEE) and Figures 2 and 7 (TTE):

- Catheter inlet area 3.5 cm below the aortic valve
- Catheter outlet area well above the aortic valve (frequently not visible on TEE or TTE images)
- Catheter angled toward the left ventricular apex away from the heart wall and not curled up or blocking the mitral valve

2—IMPELLA CATHETER TOO FAR IN LEFT VENTRICLE

If the Impella Catheter is positioned too far into the left ventricle, the patient will not receive the benefit of Impella Catheter support. Blood will enter the inlet area and exit the outlet area within the ventricle. Obstruction of the Impella Catheter inlet area can lead to increased mechanical forces on blood cell walls and subsequent hemolysis, which often presents as dark or blood-colored urine. If the Impella Catheter is too far into the left ventricle, echocardiography will likely show the following, as depicted in Figure 4 (TEE), and Figure 8 (TTE):

- Catheter inlet area more than 4 cm below the aortic valve
- Catheter outlet area across or near the aortic valve

3—IMPELLA CATHETER INLET IN THE AORTA

If the inlet area of the Impella Catheter is in the aorta, the patient will not receive the benefit of Impella Catheter support. The catheter will pull blood from the aorta rather than the left ventricle. In addition, suction is possible if the inlet area is against the wall of the aorta or valve sinus. If the inlet area of the Impella Catheter is in the aorta, echocardiography will likely show the following, as depicted in Figure 5 (TEE), Figure 9 (TTE):

- Catheter inlet area in aorta or near the aortic valve
- Catheter pigtail too close to the mitral valve

4—IMPELLA CATHETER IN PAPILLARY MUSCLE

If the inlet area of the Impella Catheter is too close to or entangled in the papillary muscle and/or subannular structures surrounding the mitral valve, it can affect mitral valve function and negatively impact catheter flow. If the inlet area of the catheter is lodged adjacent to the papillary muscle, the inflow may be obstructed—resulting in suction alarms. This positioning is also likely to place the outlet area too close to the aortic valve, which can cause outflow at the level of the aortic valve, with blood streaming back into the ventricle—resulting in turbulent flow and hemolysis. If the Impella Catheter is too close to or entangled in the papillary muscle, echocardiography will likely show the following, as depicted in Figure 6 (TEE) and Figure 10 (TTE):

- Catheter pigtail in papillary muscle
- Catheter inlet area more than 4 cm below the aortic valve or lodged between papillary muscle and the myocardial wall
- Catheter outlet area too close to the aortic valve

TRANSESOPHAGEAL AND TRANSTHORACIC ECHOCARDIOGRAPHIC IMAGES

The following images depict trans-esophageal echocardiographic (TEE) and transthoracic echocardiographic (TTE) images of four Impella Catheter positions.

Four Transesophageal Echocardiographic (TEE) images

Figure 3. Correct Impella Catheter Position (TEE)

- Catheter inlet area 3.5 cm below the aortic valve
- Catheter outlet area well above the aortic valve
- Catheter angled toward the left ventricular apex away from the heart wall and not curled up or blocking the mitral valve


Figure 4. Impella Catheter Too Far into Left Ventricle (TEE)

- Catheter inlet area more than 4 cm below the aortic valve
- Catheter outlet area across or near the aortic valve


Figure 5. Impella Catheter Inlet in Aorta (TTE)

- Catheter inlet area in aorta or near the aortic valve
- Catheter pigtail too close to the mitral valve


Figure 6. Impella Catheter in Papillary Muscle (TTE)

- Catheter pigtail in papillary muscle
- Catheter inlet area more than 4 cm below the aortic valve or lodged between papillary muscle and the myocardial wall
- Catheter outlet area too close to the aortic valve

Four Transthoracic Echocardiographic (TTE) images


Figure 7. Correct Impella Catheter Position (TTE)

- Catheter inlet area 3.5 cm below the aortic valve
- Catheter outlet area well above the aortic valve
- Catheter angled toward the left ventricular apex away from the heart wall and not curled up or blocking the mitral valve


Figure 8. Impella Catheter Too Far into Left Ventricle (TTE)

- Catheter inlet area more than 4 cm below the aortic valve
- Catheter outlet area across or near the aortic valve


Figure 9. Impella Catheter Inlet in Aorta (TTE)

- Catheter inlet area in aorta or near the aortic valve
- Catheter pigtail too close to the mitral valve


Figure 10. Impella Catheter in Papillary Muscle (TTE)

- Catheter pigtail in papillary muscle
- Catheter inlet area more than 4 cm below the aortic valve or lodged between papillary muscle and the myocardial wall
- Catheter outlet area too close to the aortic valve

COLOR DOPPLER ECHOCARDIOGRAPHY

When moving a patient supported with an Impella Catheter, it is important to monitor catheter migration. Adding color Doppler to an echo is another way to verify catheter position. If the Impella Catheter is correctly positioned, a dense mosaic pattern of turbulence will appear above the aortic valve near the outlet area of the catheter—as shown in Figure 11. If, however, the echocardiogram reveals a dense mosaic pattern of turbulence *beneath* the aortic valve—Figure 12, this likely indicates that the outlet area of the catheter is in the wrong position, that is, the catheter is too far into the ventricle or entangled in papillary muscle.

Note: If using transesophageal echocardiography [TEE], look for the mosaic patterns in the same locations relative to the aortic valve and Impella Catheter outlet area.


Figure 11. Correct Impella Catheter Position (Color Doppler TTE)


Figure 12. Incorrect Impella Catheter Position (Color Doppler TTE)

POST-INSERTION POSITIONING (PIP) CHECKLIST

Completing the steps shown in the following post-insertion positioning checklist can help to ensure proper position of the Impella Catheter following insertion. Pay particular attention to positioning after the patient is moved from the operating room or catheterization laboratory.

1. Remove slack in the Impella Catheter by increasing P-level to AUTO or P-8 (or P-9 for the Impella 5.0 or Impella LD) and align the catheter against the lesser curvature of the aorta (rather than the greater curvature).
2. Use fluoroscopy to verify that the slack has been removed.
3. Verify that the Impella Catheter inlet area is optimally positioned 3.5 cm below the aortic valve.
4. Return to previous P-level.
5. Secure the Impella Catheter at a firm external fixation point in the groin area.

UNDERSTANDING & MANAGING IMPELLA CATHETER POSITION ALARMS


The AIC continuously monitors the catheter based on the placement signal and the motor current.

- Placement signal:
Is the signal characteristic of aortic or ventricular pressure ?
(for Impella 2.5 or Impella CP)
- Is it pulsatile or flattened ? (for Impella 5.0 or Impella LD)?
- Motor current:
Is the signal "pulsatile" or "flattened"?


If the system alarms with one of the positioning alarms described in this section, echocardiography imaging is the best method for confirming position. You can also use TEE, TTE, or fluoroscopy. If the Impella Catheter is either partly, (just the pigtail) or completely, in the ventricle, reposition the catheter under imaging guidance. If the Impella Catheter is completely out of the ventricle, do not attempt to reposition the catheter across the valve without a guidewire. The following describes possible placement conditions and the associated signal characteristics and alarm messages, as well as actions to take for each. **Note:** When using Impella CP with SmartAssist or Impella 5.5 with SmartAssist with a non-optical AIC placement monitoring is not available. Monitor pump position with imaging

CORRECT POSITION

If the Impella Catheter is in the correct position, the placement screen will appear as shown in the first image below for the Impella 2.5 and Impella CP, and in the second image for the Impella 5.0 and Impella LD.

Correct Impella CP Catheter Position (similar for Impella 2.5)


Correct Impella 5.0 Catheter Position (similar for Impella LD)

IMPELLA 2.5 OR IMPELLA CP CATHETER FULLY IN VENTRICLE

If the Impella 2.5 or Impella CP Catheter is fully in the ventricle, the following alarm appears:

Impella Position In Ventricle

In this situation, the placement screen will appear as shown below.

Impella CP Catheter Fully in Ventricle (similar to Impella 2.5)

Note: Yellow Flow Rate

If the Automated Impella Controller detects incorrect or unknown catheter position, or if placement monitoring is suspended, the flow rate in the flow area appears in yellow as shown in above.

Actions to take:


1. Under fluoroscopic or echocardiographic guidance, reduce the P-level to P-2 and carefully pull back the Impella Catheter until the aortic waveform signal is showing.
2. When you see the aortic waveform signal, pull the catheter back an additional 4 cm

IMPELLA 2.5 OR IMPELLA CP CATHETER COMPLETELY IN THE AORTA OR INLET AND OUTLET AREAS IN VENTRICLE AND OPEN PRESSURE AREA IN AORTA

If the Impella 2.5 or Impella CP Catheter is completely in the aorta or if the inlet and outlet areas are in the ventricle and the open pressure area is in the aorta, the following alarm appears:

Impella Position Wrong

In this situation, the Placement Screen appears, as shown below.

Impella CP Catheter Completely in the Aorta or Inlet and Outlet Areas in Ventricle and Open Pressure Area in Aorta (similar for Impella 2.5)

Actions to take:

1. Reduce the P-level to P-2 until imaging is available.
2. Under fluoroscopic or echocardiographic guidance, determine the Impella Catheter position.
3. Return to preferred P-level and confirm positioning with imaging.

LOW NATIVE HEART PULSATILITY—IMPELLA 2.5 OR IMPELLA CP

When a patient has poor native ventricular function, the placement signal may remain pulsatile; however, the amplitude will be dampened.

In a situation of low native heart pulsatility, the Automated Impella Controller may not be able to determine the catheter position. You may see the following alarm appears:

Impella Position Unknown

Actions to take: Assess cardiac function.

IMPELLA CATHETER OUTLET AREA ON OR NEAR AORTIC VALVE

If the Impella Catheter outlet area is on or near the aortic valve, the catheter may be too deep in the ventricle.

Actions to take:


1. Assess and adjust Impella Catheter position under fluoroscopic or echocardiographic guidance.
2. If unsuccessful, reduce the P-level to P-2 and gently pull the catheter back 2 cm to see if the condition resolves.

CATHETER POSITION WRONG—IMPELLA 5.0 OR LD

If the Impella 5.0 or LD Catheter is fully in the ventricle or fully in the aorta, the following alarm will appear:

Impella Position Wrong

The Impella 5.0 or Impella LD System cannot differentiate between these two conditions. In this situation, the placement screen will appear as shown below.

Impella 5.0 Catheter Position Wrong (similar for Impella LD)

Actions to take:

1. Reduce the P-Level to P-2 until imaging is available.
2. Reposition Impella Catheter under fluoroscopic or echocardiographic guidance.
3. Return to preferred P-Level and confirm positioning with imaging.

LOW NATIVE HEART PULSATILITY—IMPELLA 5.0, IMPELLA LD

When a patient has poor native ventricular function, the placement signal may remain pulsatile; however, the amplitude will be dampened and both the minimum and maximum values will be greater than zero because the aortic valve does not open and the Impella 5.0 or Impella LD Catheter raises the aortic blood pressure above the ventricular pressure during systole.

In a situation of low native heart pulsatility, the Automated Impella Controller may not be able to determine the catheter position. You may see the following indication on the home screen:

Impella Position Unknown

Notice that the flow rate is displayed in yellow in the lower-left corner of the screen, indicating that the patient may not be getting the benefit of the displayed flow rate.

Actions to take:


1. Assess cardiac function.
2. If needed, confirm catheter position with echocardiography.

IMPELLA 5.0 OR IMPELLA LD CATHETER OUTLET AREA ON OR NEAR AORTIC VALVE

If the Impella 5.0 or Impella LD Catheter outlet area is on or near the aortic valve, the catheter may be too deep in the ventricle. The following alarm appears:

Impella Outflow Blocked

In this situation, the Impella Outflow Blocked warning appears as shown below.

Impella 5.0 Catheter Outlet Area on or near Aortic Valve (similar for Impella LD)

Actions to take:

1. Assess and adjust Impella 5.0 or Impella LD Catheter position under fluoroscopic or echocardiographic guidance, if available.
2. If fluoroscopic or echocardiographic guidance is not available, reduce the P-level to P-2 and gently pull the catheter back 2 cm to see if the condition resolves.

IMPELLA STOPPED

If the Impella Catheter has stopped suddenly:

1. Try to restart the catheter at P-8.
2. If the Impella does not restart at P-8, try to restart at P-2.
3. If the Impella does not restart or stops again, wait 1 minute and try to restart again.
4. If the Impella restarts, wean down to P-2 as the patient can tolerate. Under these circumstances, catheter function is not reliable and the Impella may stop again.
5. If the Impella does not restart, remove the Impella from the ventricle as soon as possible to avoid aortic insufficiency.

SUCTION

Suction may occur if the blood volume available for the Impella Catheter is inadequate or restricted. Suction limits the amount of support that the Impella Catheter can provide to the patient and results in a decrease in arterial pressure and cardiac output. It can damage blood cells, leading to hemolysis. It may also be an indicator of right heart failure.

SUCTION WITH IMPELLA 2.5 OR IMPELLA CP

If the Automated Impella Controller detects suction while running in **AUTO** mode, it automatically reduces motor speed to lower the flow rate to resolve the suction and displays the "Impella Flow Reduced" advisory alarm. If the suction is cleared, the controller returns the flow rate to the desired setting. If suction is still detected at the lowest motor speed, the controller displays the "Suction" alarm.

If the **Suction** or **Impella Flow Reduced** alarm occurs during Impella support, follow these recommended actions:

1. Ensure patient has adequate volume.
2. Reduce P-Level by one or two P-levels.
3. Check the Impella Catheter for correct positioning using imaging. Reposition the catheter by rotating or moving it into or out of the ventricle slightly. Either or both of these actions could help move the inlet of the Impella Catheter away from the interior ventricular wall.
4. Evaluate right ventricular function by assessing invasive hemodynamics and/or echocardiography.
5. Return P-level to pre-alarm setting.
If the Impella 2.5 or Impella CP has sudden low flows or suction at startup:
 - a) Remove the catheter from the patient and ensure that ACT is 250 seconds or above.
 - b) Closely inspect the inlet and outlet areas and remove any thrombus or other foreign materials.
 - c) Replace catheter.

SUCTION WITH IMPELLA 5.0 OR IMPELLA LD CATHETER

If the **Suction** alarm occurs during support with the Impella 5.0 or Impella LD Catheter, follow the recommended actions:

1. Reduce P-level by 1 or 2 levels to reduce the effects of suction.
2. Assess patient's fluid intake and output to confirm adequate volume status.
3. Check the Impella Catheter for correct positioning using imaging. Reposition the catheter by rotating or moving it into or out of the ventricle slightly. Either or both of these actions could help move the inlet of the Impella Catheter away from the interior ventricular wall.
4. Confirm right ventricular function by assessing CVP or right side function with echocardiography. If CVP is not an option, check the pulmonary artery diastolic pressure to assess the patient volume status.
5. Return the P-level to pre-alarm setting.

OPERATING THE IMPELLA CATHETER WITHOUT HEPARIN IN THE PURGE SOLUTION

The Impella Catheter is designed to be operated with a purge solution that contains heparin. Operation of the system without heparin in the purge solution has not been tested. In the event that a patient is intolerant to heparin, due to heparin-induced thrombocytopenia (HIT) or bleeding, physicians should use their clinical judgment to assess the risks versus benefits of operating the Impella System without heparin.

If it is in the best interest of the patient to operate the system without heparin, the dextrose solution is still required, and physicians should consider *systemic delivery* of an alternative anticoagulant. The Impella Catheter has not been tested with any alternative anticoagulants in the purge solution. Use of alternative anticoagulants may reduce the longevity or performance of the Impella catheter.

HEMOLYSIS

When blood is pumped, it is subjected to mechanical forces. Depending on the strength of the blood cells and the amount of force applied, the cells may be damaged, allowing hemoglobin to enter the plasma. Pumping forces can be generated by a variety of medical procedures including heart lung bypass, hemodialysis, or ventricular assist device (VAD) support. Patient conditions—including catheter position, pre-existing medical conditions, and small left ventricular volumes—may also play a role in patient susceptibility to hemolysis.

Hemolysis should be monitored during support. Patients who develop high levels of hemolysis may show signs of decreased hemoglobin levels, dark or blood-colored urine, and in some cases, acute renal failure. Plasma-free hemoglobin (PfHgb) is the best indicator to confirm whether a patient is exposed to an unacceptable level of hemolysis.

Management technique may differ depending on the underlying cause of hemolysis. The following table provides guidance for various circumstances.

Table Guide for Managing Hemolysis in Various Circumstances

Condition	Controller Indicators	Clinical Indicators	Management
Impella inlet area in close proximity to intraventricular wall	• “Impella Flow Reduced” or “Suction” alarms • Lower than expected flows	Imaging (see note)	• Reposition the catheter by rotating or moving the catheter into or out of the ventricle slightly. Either or both of these actions could help move the inlet of the catheter away from the intraventricular wall. • If repositioning will be delayed, reduce the P-level if tolerated by patient hemodynamics. Return to the previous P-level after repositioning. • Reassess position after flow rate has returned to desired target value.
Wrong pump position	• Position alarms with higher than expected flows • “Impella Flow Reduced” or “Suction” alarms with lower than expected flows • Pump outlet blocked alarms	Imaging (see note)	• Reposition the catheter by rotating or moving the catheter into or out of the ventricle slightly. Either or both of these actions could help move the inlet of the catheter away from the intraventricular wall. • If repositioning will be delayed, reduce the P-level to P-2. Return to the previous P-level after repositioning. • Reassess position after flow rate has returned to desired target value.
Higher than needed P-Level setting	• There may be no controller indicators • “Impella Flow Reduced” or “Suction” alarms	• Normal hemodynamics • Native recovery	• Reduce P-level until patient pressure starts to drop. • Slowly increase P-level.
Inadequate filling volume	• Position alarms • “Impella Flow Reduced” or “Suction” alarms • Lower than expected flows	• Low CVP • Low PCWP • Low AOP • High PA pressures • Right heart failure • High urine output • Increased bleeding or chest tube drainage	• Reduce the P-level if tolerated by patient hemodynamics. • Correct I and O balance. • Consider giving volume; additional volume will expand the end systolic ventricular volume. • Reduce PA pressure. • Improve right heart function.
Pre-existing patient conditions or other medical procedures	N/A	• Patient past medical history • Current procedures or treatments	
Note on imaging: All imaging technology represents the anatomy in two dimensions (2D). It is not possible to assess the interactions between the catheter and the intraventricular anatomy that occur in three dimensions (3D). Abiomed strongly recommends that the catheter be repositioned, even if the imaging view shows correct position.			

PLACEMENT SIGNAL LUMEN

(For Impella 2.5 & Impella CP)

BACKGROUND

The Impella 2.5 and Impella CP Catheters use a fluid-filled pressure lumen with an inlet at the proximal end of the motor housing and the pressure sensor located in the red Impella plug. The Automated Impella Controller software monitors both the pressure waveform characteristics and motor current to determine the placement of the Impella Catheter inlet and outlet areas relative to the aortic valve. The following table provides recommended standards for maintaining the placement signal.

Note: Restoring Placement Signal Quality


You may get a sensor or position alarm if you pinch the white flush valve to restore placement signal quality.

Recommended Standards for Maintenance of the Placement Signal for Impella 2.5 and CP Catheters

Periodic flushing of the placement signal lumen.	<i>Slight dampening</i> If you observe a dampened placement signal, pinch the white flush valve located on the red sidearm for a few seconds to restore the placement signal quality.
Note: Either of these actions may result in sensor or position alarms.	<i>Severe or lost pressure</i> 1. Close the roller clamp and disconnect the IV tubing connected to the red pressure sidearm. 2. Connect a syringe of saline to the port and squeeze the white flow valve as you draw negative pressure. 3. Continue aspiration of the port until blood is visualized in the syringe. 4. Disconnect the syringe and open the roller clamp until slow drips of saline exit the tubing. 5. Flood the open port of the red pressure sidearm and then reconnect. 6. Squeeze the white wings of the flow valve for 15 to 20 seconds to flush the pressure lumen to remove all blood from the pressure lumen.
Pressure bag inflation pressure	Maintain pressure bag inflation pressure between 300 mmHg and 350 mmHg.

FLUSH SOLUTION CHANGE-OUT PROCEDURE

1. Prime the new NaCl flush solution set-up and close the roller clamp.
2. Place the NaCl bag in a pressure bag and inflate to between 300 mmHg and 350 mmHg.
3. Close the roller clamp and disconnect the old flush solution connected at the red sidearm port.
4. Open the roller clamp on the new flush solution set-up until you get a slow drip.
5. Position the male luer connector over the female luer connector and fill to overflow, displacing any air, as shown below.
6. Connect and secure luer fittings.
7. Fully open the roller clamp and squeeze the white wings for approximately 5 to 10 seconds to complete the internal prime. This final prime should eliminate any risk of lost or dampened pressure caused by blood tracking into the pressure lumen during the pressure tubing change.


Displacing Air During Flush Solution Change Out Procedure

PRESSURE SENSOR DRIFT AND PLACEMENT SIGNAL NOT RELIABLE (FOR IMPELLA 5.0 AND IMPELLA LD)

MANUALLY ZEROING THE DIFFERENTIAL PRESSURE SENSOR

The electrical signal produced by the differential pressure sensor may drift over time. If you observe that the placement waveform has shifted up or down on the display, or the expected flow does not match the current P-level setting, zero the differential pressure sensor by performing the following steps:

1. Press the **MENU** key and select **Start Manual Zero**.
2. Select **OK** to confirm the decrease in P-level.
3. The controller displays: **Wait until the new P-level is reached** and then: **Calculation is running**.
4. Select **OK** to accept the setting when the controller displays the **Placement Signal Offset Adjust finished!**.
5. The Impella automatically resets to the previous P-level.

Placement Screen with MAPs

ZEROING THE DIFFERENTIAL PRESSURE SENSOR WHEN THE IMPELLA 5.0 CATHETER IS RUNNING

The controller software contains a data table listing the expected differential pressure for a given motor current when the motor speed is set to a specific value. To zero the differential pressure sensor while the Impella 5.0 Catheter is running, the software sets the motor speed and measures the motor current. Using the data table, the software determines what the measured differential pressure should be, then adjusts the signal from the differential pressure sensor so that it matches the expected value.

PLACEMENT SIGNAL NOT RELIABLE AND EFFECTS ON FLOW CALCULATIONS

If the pressure sensor fails, the controller can no longer calculate the flow rate. The controller displays a **Placement Signal Not Reliable** alarm. The placement screen displays a table with estimated flows and corresponding MAPs in the bottom-left corner of the display. The Home screen displays a yellow question mark over the heart icon and the text:

Placement Monitoring Suspended

To **Silence** this alarm, Go to > **MENU** and select > **SETTINGS/SERVICE**.

SUCTION DETECTION DURING SENSOR DRIFT OR SENSOR FAILURE

If sensor drift occurs or the pressure sensor fails, the controller can no longer detect suction. (For more information about suction, refer to the **Suction** discussion earlier in this manual.)

Note: The effectiveness of Impella 5.0 or Impella LD Catheter support can only be assessed by monitoring patient hemodynamics, cardiac imaging, and the Impella 5.0 Catheter motor current.

Signs of suction Include the following:

- A drop in the patient's arterial pressure
- Decreased output, if a cardiac monitor is in place
- Dampened or flat motor current waveforms

If imaging reveals that the suction is caused by the catheter inlet area in close proximity to the intraventricular wall, reposition the catheter as described in Table Guide for Managing Hemolysis in Various Circumstances.

If hemodynamic parameters, such as low aortic pressure or high pulmonary artery pressure, indicate suction caused by inadequate filling volume, reduce P-level and follow the management strategies described in the Table Guide for Managing Hemolysis in Various Circumstances.

ENABLING PURGE FLOW NOTIFICATIONS

The purge flow notification white alarms **Purge Flow Increased** and **Purge Flow Decreased**) are disabled by default.

To **Enable** these alarms:

1. Press **MENU** and scroll to **Settings/Services**. Press the selector knob.
2. Scroll to **Enable Purge Flow Change Notifications** and press the selector knob to enable these alarms.

DISABLING AUDIO ALARMS

The audio for the following alarms can be disabled:

- Placement Signal Not Reliable
- Placement Signal Lumen Blocked
- Suction
- Purge System Blocked / Purge Pressure High

To disable the audio:

1. Press **MENU** and scroll to "Settings/Services." Press the selector knob.
2. Highlight alarm and press the selector knob to disable the audio for the selected alarm.

SURGICAL MODE

Surgical Mode can be enabled to silence the **Impella Stopped** alarm that occurs when P-level is reduced to P-0. A white banner notification displays throughout the duration of Surgical Mode support.

To **Enable Surgical Mode**:

1. Press > **MENU** and scroll to > **Settings/Services**. Press the Selector knob.
2. Scroll to > **Enable Surgical Mode**, Press the Selector knob to enable it.

You can **Disable** Surgical Mode in one of two ways:

1. Increase P-level above P-0, or
2. Press **MENU**, select > **Settings/Services**, then select > **Disable Surgical Mode**.

TIMED DATA RECORDING

The Automated Impella Controller can hold up to 24 hours of real-time data. Once memory is full, the controller starts overwriting the old data. The timed data recording feature allows you to permanently save real-time operating data for later analysis. Timed data recording is automatically turned on during certain alarm conditions to capture data for analysis. You can also manually turn on the feature at any time to capture data for later analysis.

To manually access the timed data recording feature:

1. Press **MENU** and scroll to "Start Data Snapshot." Press the selector knob.
2. The controller records data for a predefined period of 10 minutes.

OPERATING THE IMPELLA CATHETER IN ELECTROMAGNETIC FIELDS

The Impella Catheter contains a permanent magnet motor that emits an electromagnetic field. This field may produce electromagnetic interference with other equipment. In addition, other equipment that emits a strong electromagnetic field may affect the operation of the Impella Catheter motor.

ELECTROANATOMIC MAPPING (EAM) SYSTEMS

Examples of EAM Systems

- CARTO® 3 System and
- CARTO® XP Navigation System (Biosense Webster, Inc.)

The electromagnetic field emitted by the Impella Catheter may produce interference with the magnetic location detection component of the electroanatomic mapping (EAM) system, particularly when the mapping catheter is close to the Impella Catheter motor. For example, mapping in the right or left ventricular outflow tracts places the mapping catheter in close proximity to the Impella Catheter motor in the ascending aorta.

Electromagnetic interference may appear as:

- Instability in the displayed location of the mapping catheter
- Magnetic interference errors generated by the electroanatomic mapping system

When operating the Impella Catheter in the presence of an EAM system, use P-level mode. Operate the Impella Catheter at P-1–P-5 or P-7.

The motor speeds at these P-levels cause the least interference. Best performance is observed when the Impella Catheter motor is at least 3 cm from the sensors in the mapping catheter. If you suspect interference, follow the troubleshooting steps in Troubleshooting When Operating the Impella Catheter in the Presence of an EAM System table.

Troubleshooting When Operating the Impella Catheter in the Presence of an EAM System

Observation	Actions
Interference with the magnetic location detection component of the EAM system	1. Check for and address other sources of interference. 2. Reposition the Impella Catheter to ensure that the Impella motor is at least 3 cm from the sensors in the mapping catheter; however do NOT pull the inlet area out of the left ventricle. 3. Ensure that the Impella Catheter is operating at P-1–P-5 or P-7, as these P-levels cause the least interference.

MAGNETIC NAVIGATION SYSTEMS (MNS)

Example of MNS

Stereotaxis Niobe® Magnetic Navigation System (Stereotaxis)

When initiating Impella Catheter support in the presence of a magnetic navigation system (MNS), follow the steps below:

1. Insert the Impella Catheter following the steps outlined in section 5 of this manual.
2. Place the MNS magnets in the "Reduced" or "Stowed" position.
3. Start the Impella Catheter in the manner described in section 5 of this manual. Increase P-level to P-5.
4. Place the MNS magnets in the "Navigate" position and proceed with magnetic navigation.

Keep operating the Impella Catheter at a P-level of at least P-5 when the MNS magnets are in the "Navigate" position. If the P-level falls below P-5, the Impella Catheter may stop running. To resume operation, follow the steps in the Troubleshooting When Operating the Impella Catheter in the Presence of a MNS System table.

During magnetic navigation of the mapping catheter, the motor current of the Impella Catheter may temporarily increase to the point that the catheter stops running. The table below explains how to resume operation. When the MNS magnets are in the "Navigate" position, the displayed Impella Catheter flow may be artificially elevated. To accurately assess the flow rate, note the displayed flow when the magnets are in the "Stowed" position.

Troubleshooting When Operating the Impella Catheter in the Presence of a MNS System

Observation	Actions
Unable to start Impella or Impella stops running	1. Place the MNS magnets in the "Reduced" position and attempt to start the Impella Catheter. 2. If the Impella Catheter does NOT start with the magnets in the "Reduced" position, place the magnets in the "Stowed" position and start the Impella Catheter. 3. Increase the Impella Catheter P-level to P-5 or higher. 4. Place the MNS magnets in the "Navigate" position and proceed with magnetic navigation.
MNS magnets: "Navigate" Displayed flow seems too high	The Impella Catheter displayed flow will be artificially elevated when the MNS magnets are in the "Navigate" position.
MNS magnets: "Stowed" Displayed flow drops	The displayed flow will be accurate when the MNS magnets are in the "Stowed" position.

TRANSFERRING FROM THE AUTOMATED IMPELLA CONTROLLER (AIC) TO A NEW AIC

TRANSFER STEPS

A backup Automated Impella Controller should be available at all times when a patient is on support. In the event that the controller fails, follow the steps below to transition the Impella Catheter to the backup controller.

1. Confirm that the backup controller is powered on and ready.
2. Press **PURGE MENU** on the original controller, select Change Purge Fluid Bag.
3. Press the **START** soft button to initiate the bolus. **Note:** the console will read "Please wait" while system completes the bolus. Do NOT change and spike purge fluid bag prior to next step.

4. Disconnect the yellow luer connector from the Impella Catheter to release the pressure in the purge cassette.
5. Transfer the purge cassette and purge solution from the original controller to the backup controller.
6. Reconnect the yellow luer connector to the Impella Catheter.
7. Remove the white connector cable from the original controller and plug it into the catheter plug on the front of the backup controller.
8. Once the Impella Catheter is connected to the backup controller, wait for a message to appear on the screen asking you to confirm re-starting the Impella Catheter at the previously set P-level.
9. Press **OK** within 10 seconds to confirm restarting the Impella Catheter at the previously set P-level.
10. If the message to restart the Impella Catheter does not appear within 30 seconds, restart the Impella Catheter using the **FLOW CONTROL** soft button.

Note: Impella CP with SmartAssist and Impella 5.5 with SmartAssist Catheters will run on a non-Optical AIC. Sensor functionality **will not work** in this case.

CHANGE PURGE FLUID TO OBTAIN ACCURATE PURGE VALUES

To get accurate purge values after changing to a backup controller, perform the Change Purge Fluid bag procedure and replace the purge fluid bag. Go to > **MENU**, select > **Change Purge Fluid**, and follow the instructions given on-screen.

PATIENT MANAGEMENT CHECKLIST FOLLOWING TRANSFER OF SUPPORT

After transferring patient support to or from the Automated Impella Controller, perform each of the following patient management checklist items:

1. Confirm Impella Catheter placement using echocardiography.
2. Tighten the Tuohy-Borst valve (tighten all the way to the right) on the Impella Catheter to prevent catheter migration.
3. For the Impella 2.5 and Impella CP Catheters, attach a saline pressure bag pressurized to 300-350 mmHg to the red sidearm and complete the "Transfer to Standard Configuration" procedure under the **PURGE MENU** option if not already completed.

EMERGENCY SHUTDOWN PROCEDURE

In the unlikely event that the Automated Impella Controller software stops responding, follow the procedure below to restart the controller.

1. Press and hold the power switch for 30 seconds.
2. An "Emergency Shutdown Imminent" alarm will sound at 15 seconds.
3. The controller will shut down after 30 seconds.
4. Restart the controller.

QUESTIONS OR CONCERNS?

Contact the local Abiomed team

or

Call the 24-hour Clinical Support Line:

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Abiomed Europe GmbH