Impella RP[®]

with the Automated Impella Controller

Circulatory Support System

**USER MANUAL
GEBRAUCHSANWEISUNG
MANUEL D'UTILISATION
MANUAL DEL USUARIO
MANUALE D'USO
ANVÄNDARHANDBOK
BRUGERVEJLEDNING
GEBRUIKERSHANDLEIDING**

Recovering hearts. Saving lives.™

TABLE OF CONTENTS

English.....	5
Deutsch.....	17
Français.....	31
Español.....	43
Italiano.....	57
Svenska.....	71
Dansk.....	83
Nederlands.....	95

INTENDED AUDIENCE

This manual is intended for use by cardiologists, surgeons, nurses, perfusionists, and cardiac catheterization lab technicians who have been trained in the use of the Impella Circulatory Support System.

WARNINGS

Warnings alert you to situations that can cause death or serious injury. The dark grey symbol  appears before warning messages.

-  The Impella RP System is intended for use only by personnel trained in accordance with the Abiomed Training Program.
-  Fluoroscopy is required to guide placement of the Impella RP Catheter. The small placement guidewire must be reliably observed at all times.
-  Be sure that the stopcock on the repositioning sheath is always kept in the closed position. Significant bleed back can result if the stopcock is open.
-  Avoid manual compression of the inlet, outlet, or sensor areas of the cannula assembly.
-  The sterile components of the Impella RP System can be used only if the sterilization indicators show that the contents have been sterilized, the packaging is not damaged, and the expiration date has not elapsed.
-  Do **NOT** resterilize or reuse the Impella RP Catheter. It is a disposable device and is intended for single use only. Reuse, reprocessing, reinserting through the introducer, or resterilization may compromise the structural integrity of the catheter and/or lead to catheter failure which, in turn, may result in patient injury, illness, or death.
-  Retrograde flow will occur from the pulmonary artery back into the inferior vena cava if the Impella RP Catheter is set at P-level P-0.
-  Do **NOT** use saline in the purge system.
-  Do **NOT** use an Impella RP System if any part of the system is damaged.
-  To prevent the risk of explosion, do **NOT** operate the Impella RP System near flammable anesthetics.
-  If at any time during the course of support with the Impella RP Catheter, the Automated Impella Controller alarms "Purge Pressure Low" or "Purge System Open," follow the instructions presented in the Automated Impella Controller manual.
-   MRI Unsafe - Do **NOT** subject a patient who has been implanted with an Impella RP Catheter to magnetic resonance imaging (MRI). The strong magnetic energy produced by an MRI machine may cause the Impella RP System components to stop working, and result in injuries to the patient. An MRI may also damage the electronics of the Impella RP System.
-  Cardiopulmonary support (CPR) should be initiated immediately per hospital protocol if indicated for any patient supported with the Impella RP Catheter. When initiating CPR, reduce the Impella RP Catheter flow rate. When cardiac function has been restored, return flow rate to the previous level and assess the placement signal on the controller.

-  During defibrillation, do **NOT** touch the Impella RP Catheter, cables, or Automated Impella Controller.
-  Avoid overinserting the Impella RP Catheter and possibly impinging the catheter tip against the walls of the vasculature, atrium, or ventricle.
-  To reduce the risk of cardiac or vascular injury (including ventricular perforation) when advancing or torquing the Impella RP Catheter, adjustments should be performed under imaging guidance.
-  Do **NOT** advance or withdraw the Impella RP Catheter against resistance without using fluoroscopy to determine the cause of the resistance. Doing so could result in separation of the catheter or guidewire tip, damage to the catheter or vessel, or perforation.
-  To reduce the possibility of fibers being drawn into the Impella, customers should avoid exposing the inlet and cannula section of the Impella Heart Pumps to any surfaces or fluid baths where the device can come into contact with loose or floating fibers.
-  To reduce the risk of cardiac or vascular injury (including perforation) when manipulating the heart during cardiac surgery, evaluate the position of the pump using imaging guidance prior to manipulating the heart, and monitor position.
-  To reduce the risk of cardiac or vascular injury (including ventricular perforation) physicians should exercise special care when inserting the Impella Catheter in patients with complex anatomy. This includes patients with known or suspected: decreased ventricular cavity size, ventricular aneurysms, thin-walled ventricles due to chronic dilation, congenital heart disease, or compromised cardiac tissue quality.
-  Physicians should exercise special care when inserting the Impella Catheter during active Cardiopulmonary Resuscitation (CPR). In addition, active CPR maneuvers may change the position of the Impella device, introducing the risk of cardiac or vascular injury (including ventricular perforation). Check that the pump is positioned correctly in the right ventricle after CPR with chest x-ray guidance.
-  Thrombus formation or deposits on indwelling central venous lines or cannulas (i.e. hemodialysis catheters, PA catheters, ECMO) may break free and enter into the Impella RP inlet, resulting in reduced flow, loss of support, or hemolysis. Assess the risk for extraluminal thrombus on indwelling lines placed prior to initiation of support.
-  Anti-coagulate patients as needed to maintain recommended ACT (160- 180s), in particular when indwelling central venous lines or cannulas (i.e. hemodialysis, PA catheters, ECMO) are present. ACT below this level may increase the risk of thrombus formation or deposition. If either internal thrombus forms within or external thrombus deposits in the Impella RP, this may result in reduced flow, loss of support, or hemolysis.
-  Do not use alcohol or alcohol-based products for cleaning.
-  Operating the pump at a high flow rate while inserting, manipulating, or removing guidewires, indwelling central venous lines or devices poses a risk of the device tip interacting with the pump inlet. This may cause damage to the interacting device, or loss of support. Reduce the P-level to P-2 and use imaging guidance during insertion, manipulation and removal of guidewires, indwelling lines or devices.

CAUTIONS

Cautions indicate situations in which equipment may malfunction, be damaged, or cease to operate. The light gray symbol  appears before caution messages.

Handle with care. The Impella RP Catheter can be damaged during removal from packaging, preparation, insertion, and removal. Do **NOT** bend, pull, or place excess pressure on the catheter or mechanical components at any time.

Patients with tricuspid or pulmonary valve stenosis or insufficiency, and patients with prosthetic tricuspid or pulmonary valves, may be compromised by the use of the Impella RP Catheter.

Use only original accessories and replacement parts supplied by Abiomed.

Do **NOT** use damaged or contaminated connector cables.

To prevent device failure, do **NOT** start the Impella RP Catheter until the placement guidewire has been removed.

Do **NOT** remove the Impella RP Catheter over the length of the placement guidewire.

When replacing the purge cassette, the replacement process must be completed within 90 seconds. The Impella RP Catheter may be damaged if replacement takes longer than 90 seconds.

To prevent malfunction of the Automated Impella Controller, avoid long-term exposure to direct sunlight and excessive heat (40°C).

To prevent overheating and improper operation, do **NOT** block the cooling vents of the Automated Impella Controller while it is operating.

Do **NOT** kink or clamp any part of the Impella RP Catheter.

Do **NOT** use the Impella RP Catheter with a damaged or kinked introducer. Replace the introducer if a kink is observed.

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella Controller will operate for at least 60 minutes after the batteries have been fully charged.

Have a backup Automated Impella Controller, purge cassette, connector cable, and Impella RP Catheter available in the unlikely event of a device failure.

Operation of the system without heparin in the purge solution has not been tested. In the event that a patient is intolerant to heparin, due to heparin-induced thrombocytopenia or bleeding, physicians should use their clinical judgment to assess the risks versus benefits of operating the Impella System without heparin. If it is in the best interest of the patient to operate the system without heparin, the glucose solution is still required, and physicians should consider systemic delivery of an alternative anticoagulant. Do **NOT** add any alternative anticoagulant (such as a direct thrombin inhibitor) to the purge fluid. The Impella Catheter has not been tested with any alternative anticoagulants in the purge solution.

Insertion through the left femoral vein may result in repeated attempts to properly position the Impella RP, which could cause excessive manipulation and pump damage. As a result, left femoral insertion should be avoided whenever possible.

Do not insert the repositioning sheath into the peel-away sheath during catheter insertion; slide the repositioning sheath towards the blue pump handle prior to inserting the device.

Never remove the guidewire when increased resistance is met. Determine the cause of resistance using fluoroscopy and verify the repositioning sheath is not inserted into the peel-away sheath.

Benefits of Impella RP in salvage patients have not been proven.

INDICATIONS

The Impella RP System (percutaneous pump for right ventricular support) is intended for clinical use in cardiology, in cardiac surgery, and intensive care unit for up to 14 days for the following indications, as well as others:

- Acute or transient reduction of the right ventricular function (eg, postcardiotomy low output syndrome)
- Cardiogenic shock as a consequence of a posterior myocardial infarction with right ventricular heart failure
- Right heart support during coronary beating heart bypass surgery, especially for patients with a reduced preoperative cardiac output or for patients having a high risk of developing a postoperative low output syndrome for other reasons
- Right ventricular heart failure after implantation of a left ventricular assist device
- Therapy unresponsive arrhythmias with a reduction of right ventricular output
- Heart failure and/or cardiogenic shock as a consequence of refractory ventricular arrhythmias, as well as a consequence of sustained supraventricular arrhythmias, causing haemodynamic compromise.

CONTRAINDICATIONS

- Arteriosclerosis, in particular calcification or other disorders of the pulmonary artery wall
- Mechanical valves, severe valvular stenosis or valvular regurgitation of the tricuspid valve or pulmonary valve
- Mural thrombus of the right atrium or vena cava
- Anatomic conditions precluding insertion of the pump
- Other illnesses or therapy requirements precluding use of the pump
- Presence of a vena caval filter or caval interruption device, unless there is clear access from the femoral vein to the right atrium that is large enough to accommodate a 22 Fr catheter

POSSIBLE COMPLICATIONS

There are risks of complications with every procedure using a blood pump. These include among others:

- Arrhythmia (atrial fibrillation, ventricular fibrillation and/or tachycardia)
- Bleeding
- Cardiac or vascular injury (including ventricular perforation)
- Cardiac tamponade
- Cardiogenic shock
- Device malfunction
- Embolism
- Endocardial injuries as a result of attachment of the pump due to suction
- Hemolysis
- Immune reaction
- Insertion site infection
- Patient dependency on the pump after use for support
- Positioning problems causing hemolysis or reduced hemodynamic support
- Pulmonary valve insufficiency
- Sepsis
- Thrombocytopenia
- Tricuspid valve injury
- Venous thrombosis.

OVERVIEW

The Impella RP Catheter is an intracardiac microaxial blood pump that supports a patient's pulmonary circulation. The Impella RP Catheter is inserted percutaneously through the femoral vein and into the pulmonary artery.

When properly positioned, the Impella RP Catheter delivers blood from the inlet area, which sits in the inferior vena cava, through the cannula, to the outlet opening in the pulmonary artery. Physicians and device operators monitor Impella RP Catheter function on the display screen of the Automated Impella Controller.

The intent of the therapy with the Impella RP System is to provide a percutaneous circulatory support system to restore normal right heart hemodynamics, reduce right ventricular work, and allow the right heart time to potentially recover adequate contractile function or to be bridged to the next therapy.

This section describes the components of the Impella RP Catheter and the Automated Impella Controller, as well as the accessory components.

REUSABLE SYSTEM COMPONENTS

The Impella RP System consists of the following reusable components:

- Automated Impella Controller — provides the user interface, alarm indications, and portable battery
- Automated Impella Controller cart — for easy transport of the Automated Impella Controller

SINGLE-USE SYSTEM COMPONENTS

The Impella RP System also includes the following single-use components:

- Impella RP Catheter
- Purge cassette
- Introducer kit
- Connector cable
- 0.027 inch, 260 cm placement guidewire

SYSTEM CONFIGURATION

Figure 1 illustrates how the Automated Impella Controller connects to the Impella RP Catheter and accessory components.


Figure 1 Automated Impella Controller, Impella RP Catheter, and Accessories

IMPELLA RP CATHETER

The Impella RP Catheter is an intracardiac microaxial blood pump that delivers up to 4.0 liters of blood per minute from the inferior vena cava into the pulmonary artery. Figure 2 illustrates the Impella RP Catheter.


Figure 2 Impella RP Catheter

DIFFERENTIAL PRESSURE SENSOR

The Impella RP Catheter has an electronic differential pressure sensor located at the proximal end of the cannula. The purpose of the pressure sensor is to generate the placement signal used to calculate the flow generated by the Impella RP Catheter.

The pressure sensor is a flexible membrane integrated into the cannula. One side of the sensor is exposed to the blood pressure on the outside of the inlet area and the other side is exposed to the pressure of the blood inside of the cannula. The sensor generates an electrical signal proportional to the difference between the pressure outside the inlet area and the pressure inside the cannula. This signal is displayed on the Automated Impella Controller as the placement signal.

PURGE CASSETTE

Do **NOT** use saline in the purge system.

The purge cassette delivers rinsing fluid to the Impella RP Catheter. The purge fluid (typically 5% glucose solution) flows from the purge cassette through the catheter to the microaxial blood pump to prevent blood from entering the motor. When the purge cassette is properly installed in the Automated Impella Controller, the Abiomed logo is upright and facing you.

DISCARD THE Y CONNECTOR

Disconnect and discard the Y connector from the purge cassette tubing. For the Impella RP System, the yellow luer on the end of the purge tubing connects directly to the yellow luer on the Impella RP Catheter.




Figure 3 Purge Cassette

ACCESSORIES

Table 1 illustrates and describes the accessories used with the Impella RP Catheter and Automated Impella Controller.

Table 1 Impella RP Catheter and Automated Impella Controller Accessories

Component	Description
	The white connector cable connects the Impella RP Catheter to the Automated Impella Controller. Clips on the cable are used to secure the purge tubing to the cable. - The socket at the black end of the cable connects to the blue Impella plug. - The white plug at the opposite end of the cable is inserted into the blue catheter plug on the front of the Automated Impella Controller.
White Connector Cable	
	The introducer kit is used to place the Impella RP Catheter. It contains: 23 Fr peel-away introducer with dilator 8 Fr, 12 Fr, 16 Fr, and 20 Fr supplemental dilators 0.035 inch x 150 cm guidewire (not shown)
Introducer kit	
	A 0.027 inch, 260 cm placement guidewire is used for the placement of the Impella RP Catheter.
0.027 inch, 260 cm Placement Guidewire	
	Hospital Provided: Glucose solution (typically 5% glucose in water with Heparin (50 IU/ml) is used as the purge fluid through the Impella RP Catheter.
Glucose Solution	
	The Automated Impella Controller cart holds the Automated Impella Controller. The cart has wheels for easy transport of the controller and a storage basket. (For more information, including assembly instructions, refer to the Automated Impella Controller cart instructions for use.)
Automated Impella Controller Cart	

IMPELLA RP PATIENT SELECTION

Benefits of Impella RP in salvage patients have not been proven.

The Impella RP patient indications and its contraindications are listed in respective sections of this Impella RP IFU document. In order to guide Impella RP users to optimal patient outcomes, Impella RP patient selection guidance is an important element. Noted below, in the following figures, are definitions and a flow chart for enabling patient selection for Impella RP based on the premarket studies that will result in optimal outcomes.

Please note in the Impella RP Best Practice Selection Algorithm that the “best practices pathway” is highlighted in dark grey and is derived from US FDA PMA studies. Patients who follow the dark grey pathway should have the best chance to benefit from Impella RP. The alternate pathway (“the light gray pathway”) represents patients who fall outside of the US FDA PMA studies guidelines, but still can be served with Impella RP. These patients are likely salvage patients and the benefit of Impella RP in these patients has not been proven.

The Exclusion Criteria “Checklist” is another useful tool for guiding optimal patient selection for Impella RP. Like the Impella RP Best Practice Selection Algorithm, this list should be used to identify patients who are most likely to benefit from Impella RP and also to identify those who may not.

RIGHT VENTRICULAR FAILURE Cardiac Index <2.2 L/min/m² despite continuous High Dose Inotropes (As defined in the RECOVER RIGHT Clinical Trial) AND ANY OF THE FOLLOWING: - CVP¹ >15 mmHg or - CVP/PCWP¹ or LAP >0.63 or - Moderate to severe global RV dysfunction on echo defined as one of the following¹: - Global RV hypokinesis or - TAPSE score of ≤14 mm or - RV diameter at base >42 mm or - RV short axis (or mid-cavity) diameter >35 mm High Dose Inotropes defined as: - Dobutamine of ≥10 µg/kg/min or equivalent for more than 15 minutes - Milrinone >120 minutes - Or administration of more than one inotrope/ vasopressor medication	Profound Cardiogenic Shock defined: - SBP <75mmHg - CI <1.3 l/min/m² - despite two or more high dose inotropes - PH <7.1 not corrected by 100 mEq NaHCO₃ - DIC - Anoxic brain Injury or CGS >24 hrs Additional Measure for Early Identification of Right Ventricular Failure $PAP/P < 1.0 \quad \left \quad PAP/P = \frac{PAS - PAD}{RA}$ Legend: - PAP/P = Pulmonary Artery Pulsatility Index - PAS = Pulmonary Artery Systolic Pressure - PAD = Pulmonary Artery Diastolic Pressure - RA = Mean Right Arterial Pressure
--	--

1. Anderson MB, et al. J Heart Lung Transplant. 2015; 34(12):1549-1560

2. Korabathina, R et al. The Pulmonary Artery Pulsatility Index Identifies Severe Right Ventricle Dysfunction in Acute Inferior Myocardial Infarction. SCAI 2012; 80:593-600

Figure 4 Best Practices - Patient Selection Guidance with Right Ventricular Failure - Impella RP Heart Pump


Figure 5 Impella RP Best Practice Selection Algorithm

Below you will find a Checklist of Exclusion Criteria for the Impella RP that can be used to optimize patient outcomes.

Treatment with the Impella RP System is appropriate for patients who develop signs of acute right ventricular failure.

The checklist below is based on the exclusion criteria for the Impella RP pre-market clinical studies. This checklist is provided to help you determine if Impella RP support is an appropriate treatment for your patient, and may predict whether your patient is likely to benefit from Impella RP support.

Pre-Insertion Assessment: Post LVAD Insertion	Pre-Insertion Assessment: Post Heart Surgery, Transplant or Acute Myocardial Infarction
☐ Time (hours) to insertion post LVAD implant _____ hrs	☐ Time (hours) to insertion post completion of Cardiac Surgery or Heart Transplant, or presentation with Acute Myocardial Infarction _____ hrs
☐ INTERMACS I patients (crash and burn with worsening lactate levels or acidosis)	☐ Patients in profound cardiogenic shock
☐ Evidence of end-organ failure (bilirubin >5 or creatinine >4 within 24 hours of implant)	☐ AMI with acute mechanical complication (VSD, ventricular rupture or pap rupture)
☐ Evidence of acute neurologic injury	☐ Unsuccessful revascularization of RCA
☐ Active infection defined as two of the following (WBC $>12,500$ or positive blood culture or fever)	☐ Active infection defined as two of the following (WBC $>12,500$ or positive blood culture or fever)
☐ RA, RV or PA thrombus	☐ RA, RV or PA thrombus
☐ Prosthetic valves in the right heart	☐ Prosthetic valves in the right heart
☐ Structural tricuspid disease	☐ Structural tricuspid disease
☐ ASD or PFO (unrepaired)	☐ ASD or PFO (unrepaired)
☐ Pulmonary valve stenosis or insufficiency	☐ Pulmonary valve stenosis or insufficiency
☐ Severe pulmonary hypertension (PAS >60 mmHg)	☐ Severe pulmonary hypertension (PAS >60 mmHg)
☐ Documented DVT and/or presence of IVC filter	☐ Documented DVT and/or presence of IVC filter
☐ Anatomic abnormalities precluding insertion	☐ Anatomic abnormalities precluding insertion
☐ PA conduit	☐ PA conduit
☐ Patients on right-sided support or ECMO	☐ Patients on right-sided support or ECMO
☐ Current Pulmonary Embolism	☐ Current Pulmonary Embolism
☐ Aortic Dissection or Marfan Syndrome	☐ Aortic Dissection or Marfan Syndrome
☐ Allergy or intolerance to contrast	☐ Allergy or intolerance to contrast
☐ HIT or sickle cell disease	☐ HIT or sickle cell disease
☐ Existing congenital heart disease that would preclude placement	☐ Existing congenital heart disease that would preclude placement
OPTIMAL CANDIDATE: Time to insertion ≤ 48 hours AND no boxes checked	OPTIMAL CANDIDATE: Time to insertion ≤ 48 hours AND no boxes checked

STARTUP

Do **NOT** use an Impella RP System if any part of the system is damaged.

The sterile components of the Impella RP System can be used only if the sterilization indicators show that the contents have been sterilized, the packaging is not damaged, and the expiration date has not elapsed.

Do **NOT** resterilize or reuse the Impella RP Catheter. It is a disposable device and is intended for single use only. Reuse, reprocessing, or resterilization may compromise the structural integrity of the catheter and/or lead to catheter failure which, in turn, may result in patient injury, illness, or death.

To prevent malfunction of the Automated Impella Controller, avoid long-term exposure to direct sunlight and excessive heat (40°C).

To prevent overheating and improper operation, do **NOT** block the cooling vents of the Automated Impella Controller while it is operating.

The Li-Ion batteries must be charged for 5 hours prior to system operation in order to meet the runtime requirement of 1 hour. Failure to do so will yield a shorter runtime. After being unplugged, the Automated Impella Controller will operate for at least 60 minutes after the batteries have been fully charged.

Have a backup Automated Impella Controller, purge cassette, connector cable, and Impella RP Catheter available in the unlikely event of a device failure.

SUPPLIES NEEDED

- Automated Impella Controller
- Impella RP Catheter and accessories
- Femoral length, balloon-tipped flow-directed catheter
- 500 cc bag of glucose solution in water for purge solution (5% recommended; 5% to 20% acceptable) with 50 IU heparin/mL

NOTE: If the purge solution is supplied in bottles, open the vent on the purge fluid spike and follow the same procedure as if supplied in bags.

CASE START

NOTE: The Impella RP Catheter is a sensitive medical device with extremely fine tolerances. In particular, the inlet and outlet areas of the catheter assembly may be damaged if subjected to strong external forces.

To reduce the possibility of fibers being drawn into the Impella, customers should avoid exposing the inlet and cannula section of the Impella Heart Pumps to any surfaces or fluid baths where the device can come into contact with loose or floating fibers.

Fluoroscopy is required to guide placement of the Impella RP Catheter. The small placement guidewire must be reliably observed at all times.

The sterile components of the Impella RP System can be used only if the sterilization indicators show that the contents have been sterilized, the packaging is not damaged, and the expiration date has not elapsed.

Avoid manual compression of the inlet, outlet, or sensor areas of the cannula assembly.

Do **NOT** remove the Impella RP Catheter over the length of the placement guidewire.

Handle with care. The Impella RP Catheter can be damaged during removal from packaging, preparation, insertion, and removal. Do **NOT** bend, pull, or place excess pressure on the catheter or mechanical components at any time.

Do **NOT** kink or clamp any part of the Impella RP Catheter.

To avoid fibers drawn into the Impella:

- Keep the Impella Heart Pump in its packaging tray until just before insertion.
- Do not attempt to run the pump in a basin of saline prior to insertion.
- Do not attempt to rinse and reinsert the device after initial insertion.
- Hold the surgical towel or 4 x 4 gauze pad away from the inflow and outflow windows, when controlling blood splatter during insertion of the Impella Heart Pump through the introducer.

CASE START

If you miss a step in the process of setting up the Impella Catheter, or if you exceed the amount of time allowed to complete a step, the Automated Impella Controller will display additional instruction screens for continuing the setup process for some steps.

For the steps discussed below, all shaded steps require sterile technique.

1. Press the **START NEW CASE** soft button from the startup screen or plug in a new Impella catheter. "Case Start" can also be selected by pressing the **MENU** soft key.
2. The controller displays the screen shown below.


Figure 6

INSERT PURGE CASSETTE

1. Open the purge cassette package onto the sterile field.
2. If included, remove and discard the Y-connector and secure the yellow luer connector on the purge tubing to the sterile field.
3. Pass the purge cassette and spike off the sterile field.
4. Spike the glucose bag/bottle.
5. Press the **NEXT** soft button to continue.
6. Open the purge cassette door by pressing the release on the left side of the controller. Insert the purge cassette into the Automated Impella Controller.
7. The purge cassette snaps into a molded compartment on the front of the controller. Follow the diagram on the inside of the purge cassette door for proper placement.
8. Slide the purge disc into the slot to the right of the purge cassette until it snaps into place. The controller will automatically begin priming the purge cassette. **NOTE: The instructions for inserting the purge disc appear if it is not snapped into place within 3 seconds of inserting the purge cassette.**
9. Extend the purge tubing and close the purge cassette door. There is sufficient room around the edges of the purge cassette door so that it will not pinch the purge tubing as it exits.
10. The controller automatically begins priming the purge cassette after it is inserted. The progress bar marks the progress of the purge cassette priming.

NOTE: Once the purge cassette is installed, be sure to close the purge cassette door to prevent the purge cassette from being dislodged accidentally.

CONNECT THE CONNECTOR CABLE

1. Remove the Impella RP Catheter from its package using sterile technique and inspect the catheter, including its connector, for damage.
2. Remove the white connector cable from its package using sterile technique.
3. Inspect the cable for damage, including damage to the connector pins at the controller end.
4. Secure the grey end of the cable to the sterile field.
5. Insert the catheter plug into the connector cable socket (grey end). The tab and the slot must be aligned during connection (see Figure 7).


Figure 7 Inserting the Catheter Plug into the Connector Cable

6. Pull back on the connection to make sure that the plug has snapped into place.
7. Snap the purge clip (located on the pressure reservoir of the clear sidearm) to the connector cable as shown in Figure 8.

NOTE: Snapping the purge clip on the pressure reservoir to the connector cable is important to prevent the tube from kinking.


Figure 8 Snapping Purge Clip to Connector Cable

8. Pass the sterile connector cable from the Impella RP Catheter off the sterile field.
9. Line up the notch on the connector cable with the notch in the blue catheter plug on the front of the Automated Impella Controller and plug the cable into the controller.
10. If you have not already done so, disconnect and discard the Y connector with the red and yellow luer from the purge tubing.
11. Connect the yellow luer on the end of the purge tubing to the yellow luer on the clear sidearm of the Impella RP Catheter as shown in Figure 9.


Figure 9 Connecting the Luer to the Impella RP Catheter

12. When the controller detects that the luer is connected, it automatically begins priming the purge lumen.

ENTER PURGE FLUID DATA

1. Enter the purge fluid information.


Figure 10

2. To select the default values displayed on the screen, press the **ACCEPT** soft button. This will select those values and automatically advance to the next screen. Note: The default purge fluid values will be the purge fluid values from the last Case.
3. To change the purge fluid information, press the **EDIT** soft buttons, scroll to the appropriate item and push the selector knob to select it or use the white arrow soft buttons. Then scroll through the values and push the selector knob or press **SELECT** to make a new selection. Press the **DONE** button to finish editing. The controller will use the default values if no other selections are made. It purge fluid values will be the purge fluid values from the last Case.
 - Purge fluid can be set to 50 mL, 100 mL, 250 mL, 500 mL, or 1000 mL.
 - Glucose concentration can be set to 5%, 10% or 20%.
 - Heparin concentration can be set to 0 IU/ml, 5 IU/ml, 6.25 IU/ml, 10 IU/ml, 12.5 IU/ml, 20 IU/ml, 25 IU/ml, 40 IU/ml, 50 IU/ml.

SECURE THE PURGE TUBING

1. To complete the setup, connect the purge tubing to the white connector cable by pushing the purge tubing into the clips attached to the white connector cable as shown in Figure 11.

Figure 11 Connecting the Purge Tubing to the Connector Cable

INSERTING THE IMPELLA RP CATHETER

NOTE: Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Each physician must evaluate the appropriateness of the procedure based on his or her medical training and experience, the type of procedure, and the type of systems used.

To reduce the possibility of fibers being drawn into the Impella, customers should avoid exposing the inlet and cannula section of the Impella Heart Pumps to any surfaces or fluid baths where the device can come into contact with loose or floating fibers.

To reduce the risk of cardiac or vascular injury (including perforation) when manipulating the heart during cardiac surgery, evaluate the position of the pump using imaging guidance prior to manipulating the heart, and monitor position.

To reduce the risk of cardiac or vascular injury (including ventricular perforation) physicians should exercise special care when inserting the Impella Catheter in patients with complex anatomy. This includes patients with known or suspected: decreased ventricular cavity size, ventricular aneurysms, thin-walled ventricles due to chronic dilation, congenital heart disease, or compromised cardiac tissue quality.

Physicians should exercise special care when inserting the Impella Catheter during active Cardiopulmonary Resuscitation (CPR). In addition, active CPR maneuvers may change the position of the Impella device, introducing the risk of cardiac or vascular injury (including ventricular perforation). Check that the pump is positioned correctly after CPR with chest x-ray guidance.

To reduce the risk of cardiac or vascular injury (including ventricular perforation) when advancing or torquing the Impella RP Catheter, adjustments should be performed under imaging guidance.

Thrombus formation or deposits on indwelling central venous lines or cannulas (i.e. hemodialysis catheters, PA catheters, ECMO) may break free and enter into the Impella RP inlet, resulting in reduced flow, loss of support, or hemolysis. Assess the risk for extraluminal thrombus on indwelling lines placed prior to initiation of support.

Anti-coagulate patients as needed to maintain recommended ACT (160- 180s), in particular when indwelling central venous lines or cannulas (i.e. hemodialysis, PA catheters, ECMO) are present. ACT below this level may increase the risk of thrombus formation or deposition. If either internal thrombus forms within or external thrombus deposits in the Impella RP, this may result in reduced flow, loss of support, or hemolysis.

Fluoroscopy is required to guide placement of the Impella RP Catheter. The small placement guidewire must be reliably observed at all times.

Avoid manual compression of the inlet, outlet, or sensor areas of the cannula assembly.

Do **NOT** kink or clamp any part of the Impella RP Catheter.

Handle with care. The Impella RP Catheter can be damaged during removal from packaging, preparation, insertion, and removal. Do **NOT** bend, pull, or place excess pressure on the catheter or mechanical components at any time.

IMPELLA RP INSERTION STEP-BY-STEP (23 FR PEEL-AWAY)

NOTE: All shaded steps below require sterile technique.

1. Confirm purge fluid is exiting the Impella Catheter.
2. Obtain access to the femoral vein.
3. Insert a 5–8 Fr introducer over the 0.035 inch guidewire (provided) to pre-dilate the vessel.
4. Remove the 5–8 Fr introducer over the 0.035 inch guidewire. Insert the 8 Fr, 12 Fr, 16 Fr, and 20 Fr dilators sequentially, as needed. Remove the 20 Fr dilator and insert the 23 Fr introducer with dilator. While inserting the 23 Fr introducer, hold the shaft of the introducer to advance it into the vein.
5. Administer heparin. When ACT is at least 250 seconds, remove the 23 Fr dilator.
6. Insert a flow-directed balloon-tipped catheter into the 23 Fr introducer and advance it over a guidewire into the left (preferred) or right pulmonary artery.
7. Remove the 0.035 inch guidewire, leaving the flow-directed balloon-tipped catheter in the pulmonary artery. Form a curve or bend on the 0.027 inch, 260 cm placement guidewire and then insert it.
8. Advance the placement guidewire deep into the LPA until wire prolapses.
9. Remove the flow-directed balloon-tipped catheter.
10. Wet the cannula with sterile water and backload the catheter onto the placement guidewire. One or two people can load the catheter on the guidewire.

One-person technique

- a. Advance the guidewire into the Impella Catheter and stabilize the cannula between the fingers. This prevents pinching of the outlet area. The guidewire must exit the inlet area on the inner radius of the cannula and align with the straight black line on the catheter. The cannula can be hyperextended as necessary to ensure the guidewire exits on the inner radius of the cannula.

Two-person technique

- b. The scrub assistant can help stabilize the catheter by holding the catheter proximal to the motor. This will allow the implanting physician to visualize the inner radius. The guidewire must exit the inlet area on the inner radius of the cannula and align with the straight black line on the catheter. The physician can focus on advancing the guidewire and, if the cannula needs to be hyperextended, the scrub assistant is available to assist.
11. Advance the catheter through the hemostatic valve into the femoral vein and along the placement guidewire using a fixed-wire technique. Follow the catheter under fluoroscopy, and rotate the catheter as it enters the right ventricle to direct the cannula tip upward and across the pulmonary valve. Position the outlet area of the cannula approximately 4 cm past the pulmonary valve annulus. **NOTE: While the entire pump is in the abdominal IVC, calibrate the sensor by pressing the ZERO SENSOR soft button.**
12. Remove the placement guidewire.
13. Confirm position with fluoroscopy.

USE FLUOROSCOPY FOR PLACEMENT

Impella RP Catheter performance will be compromised if correct placement cannot be confirmed. While other imaging techniques, such as transesophageal echocardiography (TEE), can help confirm the position of the Impella RP Catheter after placement, TEE does not allow visualization of the entire catheter assembly and is inadequate for reliably placing the Impella RP Catheter.

POSITIONING AND STARTING THE IMPELLA RP CATHETER

Retrograde flow will occur from the pulmonary artery back into the inferior vena cava if the Impella RP Catheter is set at P-level P-0.

To reduce the risk of cardiac or vascular injury (including ventricular perforation) when advancing or torquing the Impella RP Catheter, adjustments should be performed under imaging guidance.

When the Impella RP is properly positioned across the pulmonary valve, but is not yet running, the placement signal will be similar to a pulmonary arterial waveform. After starting the Impella RP, the amplitude of the placement signal will increase by a factor of 2 to 2.5, depending on the selected P-level.

When the Impella Catheter is not correctly placed, there is no effective unloading of the ventricle. The patient may not be benefiting from the flow rate shown on the controller.

1. Press the **START IMPELLA** soft button.
2. Turn the selector knob to increase the P-level from P-0 to P-2.
3. Press the selector knob to select the new P-level.
4. Press the **YES** soft button to confirm after P-level is selected.
5. The catheter operation icon in the lower left corner of the screen begins rotating when the Impella RP Catheter begins to operate.
6. Increase P-level to P-9 to confirm correct and stable placement. Evaluate the catheter position and remove any excess slack. The catheter inlet area should be in the inferior vena cava and the outlet area in the pulmonary artery. Verify placement with fluoroscopy.

P-LEVELS

You can select one of ten P-levels (P-0 to P-9) as shown in Table 2. Flow rate is increased by approximately 10% with every additional P-level, but depends on preload and afterload and can vary due to suction or incorrect positioning. Select the lowest P-level that will enable you to achieve the flow rate necessary for patient support.

Table 2 P-Level Flow Rates

P -LEVEL	FLOW RATE (L/min)
P-0	0.0
P-1	0.0 – 1.2
P-2	0.0 – 1.6
P-3	0.0 – 2.0
P-4	1.3 – 2.9
P-5	1.6 – 3.1
P-6	2.4 – 3.5
P-7	3.0 – 4.0
P-8	3.4 – 4.2
P-9	3.9 – 4.4

**Flow rate depends on preload and afterload and can vary due to suction or incorrect positioning.*

Retrograde flow will occur from the pulmonary artery back into the inferior vena cava if the Impella RP Catheter is set at P-0. Retrograde flow may also occur at P-1.

At P-levels between P-1 and P-6, the Impella RP operates with a regularly recurring rapid speed pulse. This minimizes stasis and reduces the risk of thrombosis in the motor area.

SUCTION

If suction is an issue, the flow displayed on the controller may be higher than the actual Impella RP flow rate. If the suction alarm appears on the controller when the Impella RP is running at P-levels between P-7 and P-9, decrease P-level to P-6, or to P-5 or P-4, as needed, to resolve suction. If the suction alarm continues when the P-level is between P-4 and P-6, momentarily stop the Impella RP to resolve the suction issue and then restart it immediately.

USE OF THE REPOSITIONING SHEATH AND THE 23 FR PEEL-AWAY INTRODUCER

To reduce the risk of cardiac or vascular injury (including ventricular perforation) when advancing or torquing the Impella RP Catheter, adjustments should be performed under imaging guidance.

1. Flush the sidearm of the repositioning sheath located on the catheter shaft.
2. Attach a stopcock and flush the repositioning sheath prior to advancing the sheath.
3. Apply manual pressure above the puncture site and remove the 23 Fr peel-away introducer completely from the vein over the catheter shaft.
4. Grasp the two wings and bend back until the valve assembly comes apart. To do this, first stretch then snap the flexible valve mechanism that temporarily holds the two wings together. Continue to peel the two wings until the introducer is completely separated from the catheter shaft. (Note: Do NOT peel the 23 Fr peel-away introducer over the tip of the repositioning sheath.)
5. Place two dead-end caps on the repositioning sheath stopcock to prevent further usage. The sideport should not be used to give medication or draw blood because the blood could potentially clot. Pressure bags should not be connected to the sideport of the repositioning sheath. If a pressure bag is connected, the sideport must have an infusion pump or flow limiting valve in place to control the amount of fluid administered to the patient.
6. Slide the repositioning sheath over the catheter shaft and advance it into the femoral vein to the yellow eyelet.
7. Make sure there is no bleeding at the transition from the repositioning sheath to the femoral vein. Close and dress the wound.
8. Secure the repositioning sheath by suturing it to the skin using the yellow eyelet on the sheath.
9. Attach the anticontamination sleeve to the yellow section of the repositioning sheath. Lock the anchoring ring in place by turning it clockwise. Secure the catheter shaft in place by tightening the connected anchoring ring.
10. Carefully extend the anticontamination sleeve to maximum length and secure the end closest to the blue Impella plug by tightening the anchoring ring.
11. Reposition the catheter as necessary.

To prevent contamination and subsequent infections, always use sterile technique on the insertion site. Follow institutional protocols for prophylaxis of infection for patients on ventricular assist devices and indwelling lines, as well as protocols for surveillance of such patients. If device-related infection occurs, consider each patient's clinical circumstances when deciding whether to continue Impella support.

INSERTION, MANIPULATION, OR REMOVAL OF CONCOMITANT DEVICES

Guidewires, indwelling central venous lines, or devices have the potential to interact with the Impella RP inlet area during use. The following guidance is recommended:

- When inserting, manipulating, or removing guidewires, indwelling central venous lines, or devices, reduce the P-level to P-2 and use imaging guidance to limit the interaction between the distal tip of the device and the inlet area of the Impella RP.
- When manipulating or removing existing indwelling central venous lines or cannulas, thrombus formation or deposits on indwelling central venous lines or cannulas may break free. Reduce the P-level to P-2 during manipulation or removal to reduce the affinity of potential thrombus from entering the inlet area of the Impella RP.
- After insertion, manipulation, or removal is complete, return to previous P-level.

If the patient is maintained at a low flow rate (<1.5L/min) for a prolonged period of time, increase ACT to at least 250 seconds.

PATIENT WEANING

Weaning the patient from the Impella RP Catheter is at the discretion of the physician. Weaning may occur when right ventricular recovery is suspected and/or the patient is approaching the maximum duration of use for the Impella RP Catheter. It should be initiated in a step-wise manner, such as described below.

As right-side support is slowly weaned, right ventricular recovery is indicated by preservation of the normal range of left-sided cardiac output as well as by a lack of severe elevation in CVP.

The following weaning protocol is provided as guidance only.

1. Initiate the weaning process by temporarily reducing the Impella RP Catheter flow to about 2 L/min.
2. Assess right ventricular function. Small changes in right ventricular systolic function as measured by echocardiography may be accompanied by significant improvement in right side forward flow; therefore, it is important to evaluate both echocardiographic evidence of improvement as well as CVP, flow rate, and overall perfusion.
3. Record available information regarding flow rate, CVP, echo parameters, and systemic hemodynamics.
4. After 15–20 minutes at the reduced flow rate, if there are signs of right ventricular recovery and no adverse effects from reduction in flow rate, continue the weaning process by reducing flow rate as tolerated to 0.5 L/min (P1). At this flow rate there will no longer be any forward flow across the right heart.
5. If the patient is maintained at a low flow rate (<1.5 L/min) for a prolonged period, increase ACT to at least 250 seconds.

REMOVING THE IMPELLA RP CATHETER

1. Wean the patient by following the steps in the previous section.
2. Leave the Impella RP in the pulmonary artery at P-2 until ACT drops below 150
OR
Reduce the P-level to P-1, pull the catheter into the inferior vena cava (approximately 30 to 40 cm), and wait until ACT drops below 150.
3. When ACT is below 150 seconds and patient hemodynamics remain stable, decrease P-level to P-1, pull the catheter into the inferior vena cava if it is not already there, and stop the motor by reducing the P-level to P-0.
4. Remove initial mattress suture and place new mattress suture, but do not tie off.
5. Remove the Impella RP Catheter.
6. Tie off mattress suture. Apply pressure until hemostasis is achieved.
7. Disconnect the connector cable from the Automated Impella Controller and turn the controller off by pressing the power switch on the side of the controller for 3 seconds.

OPERATING THE IMPELLA CATHETER WITHOUT HEPARIN IN THE PURGE SOLUTION

The Impella Catheter is designed to be operated with a purge solution that contains heparin. Operation of the system without heparin in the purge solution has not been tested. In the event that a patient is intolerant to heparin, due to heparin-induced thrombocytopenia (HIT) or bleeding, physicians should use their clinical judgment to assess the risks versus benefits of operating the Impella System without heparin.

If it is in the best interest of the patient to operate the system without heparin, the glucose solution is still required, and physicians should consider systemic delivery of an alternative anticoagulant. **DO NOT** add any alternative anticoagulant (such as a direct thrombin inhibitor) to the purge fluid. The Impella Catheter has not been tested with any alternative anticoagulants in the purge solution.

GENERAL PATIENT CARE CONSIDERATIONS

- Do not raise the head of the bed to higher than a 30-degree angle.
- Use knee immobilizer as needed to maintain access site straight.
- Perform dressing changes per hospital protocol, using aseptic technique.
- Assess access site for bleeding and hematoma.
- Assess access site for infection.
- Be careful not to pull on the Impella RP Catheter when transferring a patient from one bed to another.
- Use care when moving or turning a patient; the Impella RP Catheter may move out of position and cause a positioning alarm.
- Assess lower extremity for DVT. The institutions should follow their own protocols on DVT prevention such as using compression devices on lower extremities and/or regular ultrasound of the legs while an indwelling vein catheter is in place.

ALARMS OVERVIEW

The Automated Impella Controller monitors various functions to determine whether specific operational parameters are within expected limits. When a parameter goes outside of its specified limits, the Automated Impella Controller sounds an alarm tone and displays an alarm message that can be viewed on the display screen on the front of the controller. The alarm tone indicates the severity of the alarm. The alarm message on the display screen is color-coded for severity and provides details on the cause of the alarm and how to resolve the alarm. After muting an alarm, if another alarm occurs it will only be heard and displayed if it is a higher priority alarm than the one that was muted.

For a full list of all Impella alarms, refer to the Automated Impella Controller Instructions for Use document.

ALARM LEVELS

Alarms are divided into three levels of severity:

- Advisory (white)
- Serious (yellow)
- Critical (red)

Table 3 Alarm Levels

Category	Description	Audible Indicator*	Visual Indicator
Advisory	Notification	1 beep every 5 minutes	Alarm header on white background
Serious	Abnormal situation. Prompt action needed.	3 beeps every 15 seconds	Alarm header on yellow background
Critical	High priority. Immediate action needed.	10 beeps every 6.7 seconds	Alarm header on red background

* Sound pressure of audible alarm indicators is >80 dBA

For some alarms, there is a short delay between the triggered event and the audible annunciation and visual display of the alarm.

MUTE ALARM FUNCTION

Pressing the **MUTE ALARM** button on the upper right of the Automated Impella Controller display screen will silence the audible alarm indicator for 2 minutes (for red or yellow alarms) or 5 minutes (for white advisory alarms). When an alarm is silenced, the words "MUTE ALARM" next to the button are replaced by the mute alarm indicator, a crossed-out bell icon.

The audible indicator will shut off if an alarm condition is resolved before you press **MUTE ALARM**. The visual message, however, will continue to be displayed, with the alarm header on a gray background, for 20 minutes or until you press **MUTE ALARM**. This allows you to identify the alarm that occurred.

TERMINOLOGY, ABBREVIATIONS, AND SYMBOLS

TERMINOLOGY AND ABBREVIATIONS

Catheter serial number	Identification number of the Impella RP Catheter; stated on the package label, on the blue Impella plug, and the Automated Impella Controller display screen
Dextrose and Glucose	The terms "dextrose" and "glucose" are used interchangeably to refer to the solution used as purge fluid for the Impella RP System
Hz	Hertz
Motor housing	Enclosure of the Impella RP Catheter motor
Pump	Central delivery unit of the Impella RP Catheter, consisting of the motor, motor housing, cannula with inlet and outlet, and pigtail at the tip
Purge Pressure	Pressure present in the Impella RP Catheter and in the infusion line
Purge System	Impella purge cassette used for rinsing the Impella RP Catheter
Retrograde Flow	Reverse flow through the cannula when the Impella RP Catheter is at a standstill (e.g., regurgitation)
V	Volt
VA	Volt ampere (Watt)

SYMBOLS

			Caution; follow instructions for use
			Defibrillator-proof type CF equipment
			Keep dry
			Storage temperature (eg, 10°C to 25°C)
			Declares conformity with Directive 93/42/EEC for medical devices or with Regulation (EU) 2017/745 on Medical Devices, and Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment
			2019-11-06 Date of manufacture (e.g., 2019-11-06)
			Protect from sunlight
			Symbol for lot designation; the manufacturer's lot designation must be stated after the LOT symbol

	123456	Abiomed part number (eg, part number 123456)
	123456	Manufacturer's serial number (eg, serial number 123456)
	2019-11-06	Use-by date (e.g., 2019-11-06)
	Do not reuse	
	Sterilized using ethylene oxide	
	Medical Device	
	Do not use if package is damaged	
	Single sterile barrier system with protective packaging inside	
	Do Not Flush	Do NOT flush with syringe
	Use glucose in the purge fluid	
	Do Not Clean with Alcohol	Do not use alcohol or alcohol-based products for cleaning.
	MR Unsafe	
	Manufacturer	
	Authorised Representative in the European Community	
	Importer	

WHITE CONNECTOR CABLE

Length	2.5 m
Service Life	Single use only

IMPELLA RP CATHETER PARAMETERS

NOTE: The Automated Impella Controller and Impella RP Catheter, including all accessories, are latex free.

Speed range	0 to 33,000 rpm
Power consumption	Less than 23 W
Voltage	Max. 20 V DC
Flow-Maximum	4.0 L/min
Purging the Impella RP Catheter	
Recommended purge fluid	5% glucose solution with heparin (50 IU per ml)
Glucose concentration	5% to 20%
Purge pressure	300 to 1100 mmHg
Purge flow	2 to 30 mL/hr
Maximum duration of use	
EU / Canada	Up to 14 days
Dimensions of Impella RP Catheter	
Length of invasive portion (without catheter)	Approx 248 mm
Diameter	Max. 7.6 mm (nom. 7.3 mm)
Classification per DIN EN 60601-1	
	Protection class I, degree of protection: CF (Automated Impella Controller and Impella RP Catheter)
Classification per directive 93/42/EEC	
	Class III
Latex free	
	Yes

RETURNING AN IMPELLA RP CATHETER TO ABIOMED

To return an Impella RP Catheter to Abiomed, contact your local Clinical Consultant for an Abiomed-approved return kit. The kit includes instructions for returning the Impella RP Catheter to Abiomed.

IMPELLA CATHETER DIMENSIONS


Figure 12

DISPOSING OF THE IMPELLA CATHETER AND ACCESSORIES

The Impella Catheter and connector cable are disposable items that must be disposed of in accordance with hospital regulations for blood contaminated materials.

The Automated Impella Controller is marked according to Directive 2012/19 EU. Devices sold within the EU can be returned to Abiomed Europe GmbH for correct disposal.

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Abiomed, Inc.